

Securos Surgical®

Instrument care and wellness product catalog



Instrument care and wellness products

Pre-cleaning

Instrument prep-clean foaming spray

Item #117470 (22 fluid ounce spray bottle)

Instrument prep-clean foaming spray is a safe and ideal method for keeping soiled instruments and scopes moist, helping to prevent bio-burden from adhering to the surface by breaking down blood, fat, protein, and carbohydrates before manual cleaning.

- Helps to prevent dried on bio-burden
- Begins breaking down blood, fat, protein, and carbohydrates
- Safe for all instruments and scopes – extends usable lifetime
- Neutral pH (non-abrasive) will not harm any metals, plastics, rubber, or corrugated tubing



Manual and ultrasonic cleaners

Power-ease instrument cleaner

Item #117464 (1 gallon)

Premium non-enzymatic detergent designed for use in all washers, washer sterilizers, ultrasonics, and for manual cleaning.

- Removes protein, blood, tissue, fat, and any organic soils by automatic washers, ultrasonic cleaner, or manual cleaning
- For use on all surgical instruments, plastics, glassware, and rubber tubing
- Fast wetting, odorless, non-sudsing, maintains high detergency and exceptional cleaning in hard, soft, cold, or hot water
- Biodegradable; contains no phosphates or NTA



Multi-enzymatic instrument cleaner

Item #117468 (1 gallon)

Best-in-class detergent with superior soil penetration and suspension; rapidly breaks down tough-to-clean medical soils including the multi-layers of bio-burden and prevents redeposition.

- Enhanced enzymatic activity combined with super detergency rapidly removes protein-rich medical soils and bio-burden
- Instantly active, constantly active, begins cleaning on contact
- Highly effective in both hard and soft water and in all water temperatures
- Neutral pH range, non-abrasive, low-foaming, free-rinsing, non-toxic, non-corrosive, and environmentally friendly



Lubricating and refinishing products

Premix and concentrate instrument lubricants

Premix item #117465 (16.9 fluid ounce spray bottle)

Premix item #117461 (1 gallon)

Concentrate item #122788 (1 gallon)

Instrument lubricant provides permanent lubrication and protection for all surgical instruments. Premix is available in ready-to-use spray and immersion applications. Instrument lubricant concentrate is available for immersion application only and needs to be diluted.

- High-grade mineral oil lubricant and rust inhibitor
- Prevents rust, spots, stains, and corrosion thus prolonging the life of the instrument
- Prevents damage due to the misuse of stiff or grinding instruments thus reducing repair and replacement costs
- Steam penetrable, steam sterilizable, and ETO sterilizable
- Non-oily, non-sticky, non-toxic, and silicone free



Instrument revitalizer

Item #117473 (16.9 fluid ounce bottle)

Instrument revitalizer is a safe, efficient, and quick revitalizing solution for stainless steel surgical instruments, trays, basins, and case carts. The product removes rust, stains, spotting, hard water scale and mineral deposits frequently encountered from sterilization. It also helps to loosen stiff joints and locks.

- Safe, efficient, and quick revitalizing solution that removes rust, stains, spotting, hard water scale and mineral deposits from stainless steel and surgical instruments
- Excellent for use in ultrasonic machines
- Routine use will enhance the life and efficiency of surgical instruments helping to loosen stiff joints and locks
- Free rinsing, non-aerosol, environmentally safe



Foaming autoclave cleaner

Item #117471 (22 fluid ounce spray bottle)

Foaming autoclave cleaner is a safe, efficient, and quick foaming revitalizing solution. It will remove all rust, spots, stains, and corrosion, as well as scale, alkaline, and mineral deposits on the inside walls of the autoclave and from stainless steel carts and trays.

- Removes rust, stains, corrosion, scale, alkaline, and mineral deposits
- Foaming spray remains on walls longer for better penetration
- Eliminates all noxious odors
- Restores original finish to autoclaves, carts, and trays



Learn more about our cleaning products, instrumentation, and available services at [securus.com](https://www.securus.com)

Orders: 800.824.3703 Service support: 877.266.3349 option #2

Premium **quality**,
uncompromising **value**

securus.com

Revised August 2025