

The DrillCover Hex System

The DrillCover Hex System is designed for small and large animal veterinarians. The DrillCover transforms a low-cost power drill into a 100% sealed and sterile surgical drilling solution, letting you perform worry-free orthopaedic surgery.

Features and Benefits

1. A fully-sealed Level 4 barrier that provides the highest protection against fluids and pathogens.
2. Paired with a robust DeWalt power drill that matches the speed, torque and feel of a surgical drill.

Cover Specifications

Barrier Level	Level 4 ANSI/AAMI PB70 barrier performance against fluid penetration
Lifetime	Linen: certified to 75 use cycles as Level 4 barrier Chuck: tested over 300 use cycles, likely significantly longer lifetime
Registration	Securos is a registered establishment with the US FDA
Sterilization	Steam autoclave 30 minutes exposure time at 121 C

Drill Specifications

Voltage	12 V battery, rechargeable on 110-120 V AC
Speed	Variable drill speed 0-1050 RPM
Torque	Drill stall-torque of 27.3 in-lbs (3.1 N-m); comparable with surgical drills
Battery	2x lithium ion, rechargeable, at least 30 minute continuous unloaded run time on full charge



"I have used the DrillCover System from Arbutus Medical and have found it to be a useful solution for practitioners who do some orthopedic surgery but cannot justify the more expensive surgical drills. The system is usable and practical for simple procedures requiring drilling or driving short pins."

- **Walter Renberg (DMV, MS, DACVS)**
Kansas State University





Endowed sponsor of the Veterinary Orthopedic Society



Orders: 800.824.3703, Service Support: 877.266.3349



securos.com



Securos Surgical®
AmerisourceBergen

How It Works

"I have used a few other drill covers in the past, and the Arbutus model is my hands-down favorite. Other drills have the chuck adapter simply passing through an un-sealed hole in the cover - not ideal from a sterility perspective. The Arbutus cover has a sealed adapter incorporated into the cover, greatly improving sterility. The cover itself is made of a material that allows very good tactile feel of the drill beneath, and the seal across the bottom (where the drill slides in) is the most secure of any that I have used. I look forward to years of using this drill for my veterinary orthopedic procedures, and would highly recommend it to anyone who needs a reasonably priced, sterile, and effective surgical drill cover."

- Dr. Ken Gummesson

Aberdeen Veterinary Hospital, Kamloops, BC

Item #	Description	Mfg #
062723	DRILL KIT COMPLETE ARBUTUS	KVI-001
062794	DRILL LINEN ARBUTUS	RP-001
062795	ADAPTOR CHUCK ARBUTUS	RP-002
062796	CHUCK KEY ARBUTUS	CHK-002

Arbutus Medical

#560-828 West 10th Ave
Vancouver, BC V5Z1L8



Arbutus Customer Service

✉ veterinary@arbutusmedical.ca



+1-778-296-2577 Option #4

Top Photo: Dr. Amy Murphy (Peace River Veterinary Clinic, Alberta)

Bottom Photo: Dr. Geoffrey Hutchinson DVM, MS, Diplomate ACVS
(Boundary Bay Veterinary Specialty Hospital, BC, Canada)



MWI® Animal Health
AmerisourceBergen



800.824.3703



mwianimalhealth.com



AmerisourceBergen®