


# Premium Sutures





We're Securos Surgical,  
the standard for  
customer experience.

## Who is Securos Surgical?

Securos is an innovator of practical solutions for veterinary surgeons. We specialize in orthopedics, but further equip the surgical suite by offering quality daily-use items including surgical instruments and suture.

## What makes us different? Why should you do business with us?

You can trust us for innovation, quality and value. And, as important as the products we bring to market, we deliver these products with the industry's most personalized service.

- Customer-focused in every aspect of what we do and every decision we make
- Problem-solving and inviting of collaboration with our customers to solve those problems
- Experts in our domain, but humble enough to constantly seek feedback
- Clear in our commitment to practical innovation, quality, value and service
- Rooted in the day-to-day reality of veterinary surgeons and general practitioners

### Orthopedics

Developed from the insights of practicing veterinary surgeons.

### Instruments

Performance you can feel the moment each instrument reaches your hand.

### Sutures

Premium quality at price points that drive your success.

## TABLE OF CONTENTS

<b>Suture Profiles</b>	<b>4-5</b>
<b>Packaging Design &amp; Comparison</b>	<b>6-7</b>
<b>Suture Needles</b>	<b>8-10</b>
Securos 300™ Needles	
Needle Types	
Cobra Black® Needles	
<b>Absorbable Suture</b>	<b>11-20</b>
Securosorb™ Quick (PGA)	
Securosorb™ (PGLA)	
Securocryl™ (PGCL)	
Securdox™ (PDO)	
Chromic Gut	
<b>Antibacterial “PLUS” Absorbable Suture</b>	<b>21-24</b>
Securosorb™ PLUS (PGA)	
Securdox™ PLUS (PDO)	
Securocryl™ PLUS (PGCL)	
<b>Non-Absorbable Suture</b>	<b>25-31</b>
Securolon™ (NYLON)	
Securopro™ (POLYPROPYLENE)	
Visipro (POLYPROPYLENE)	
SecuroSilk™ (SILK)	
<b>Suture Cassettes</b>	<b>32-33</b>
Securosorb™ (PGA)	
Securocryl™ (PGCL)	
Securdox™ (PDO)	
Chromic Gut	
Visipro (POLYPROPYLENE)	
Securomid™ (NYLON)	
<b>Comparison and Cross Reference</b>	<b>34-39</b>
Material Comparison	
SKU Quick Reference Charts	

## SUTURE PROFILES


### ABSORBABLE SUTURE



#### SECUIROSORB™ QUICK

Fast-absorbing, Coated, Braided, Polyglycolic Acid (PGA), Synthetic Absorbable Suture


**Undyed** Shelf-life: 5 years



#### SECUIRODOX™

Monofilament, Polydioxanone (PDO) Synthetic Absorbable Suture

**Violet** Shelf-life: 5 years



#### SECUIROSORB™

Coated, Braided, Polyglactin 910 (PGLA), Synthetic Absorbable Suture


**Violet or Undyed** Shelf-life: 5 years



#### CHROMIC GUT

Chromic Treated, Natural Absorbable Suture

**Beige** Shelf-life: 5 years



#### SECUIROCRYL™

Monofilament, Poly(glycolide-co-caprolactone) (PGCL) Synthetic Absorbable Suture

**Violet** Shelf-life: 4 years

### ANTIBACTERIAL ABSORBABLE SUTURE



#### SECUIROSORB™ PLUS ANTIBACTERIAL

Coated, Braided, Polyglactin 910 (PGLA), Synthetic Absorbable Suture

**Violet** Shelf-life: 3 years



#### SECUIRODOX™ PLUS ANTIBACTERIAL

Monofilament, Polydioxanone (PDO) Synthetic Absorbable Suture

**Violet** Shelf-life: 3 years




#### SECUIROCRYL™ PLUS ANTIBACTERIAL

Monofilament, Poly(glycolide-co-caprolactone) (PGCL) Synthetic Absorbable Suture

**Violet** Shelf-life: 3 years

## TENSILE STRENGTH COMPARISON CHART



# SUTURE PROFILES

## NON-ABSORBABLE SUTURE



### SECUROLON™

Monofilament Nylon  
Non-Absorbable Suture

**Black or Blue** Shelf-life: 5 years



### SECUROSILK™

Natural Silk  
Non-Absorbable Suture

**Black** Shelf-life: 5 years



### SECUROPRO™

Monofilament Polypropylene  
Non-Absorbable Suture

**Blue** Shelf-life: 5 years



### SECUROMID™

Braided, Cable, Nylon  
Non-Absorbable Suture

**White** Shelf-life: 5 years



### VISIPRO\*

Fluorescent  
Monofilament Polypropylene  
Non-Absorbable Suture

**Pink** Shelf-life: 5 years

\* VISIPRO is a trademark of Riverpoint Medical



## PACKAGING DESIGN - FEATURES & BENEFITS



### ELEGANT SUTURE BOX

- Color coded for easy identification.
- Easy open box over-wrap pull tab offers swift and painless opening.
- One-dozen size is optimal for storage after use.

### COLOR IDENTIFICATION POUCH PULL SLOT

### CONVENIENT SUTURE PEEL POUCH

- Quality Every Time. Foil pouch offers complete moisture and UV barrier.
- Easy to open Securos peel design.
- Single foil pouch unit. Reduces waste and steps to get product onto sterile field.



### PROTECTIVE FOIL POUCH EASY PEEL DESIGN

### PREMIUM SUTURE RACETRACK

- High Performance Racetrack card significantly reduces suture memory.
- Optimal Handling. Racetrack alleviates the hassle involved with tugging and straightening the suture.
- Needle Quality. Premium needle tip protected in package by foam.



### SHARPER NEEDLES LESS MEMORY OVAL RACETRACK CARD

# PACKAGING COMPARISON

Securos' packaging is specifically designed to make it a quick and smooth conversion from the leading brand to our product.




LEADING BRAND





SECUROS

PRODUCT CODE  
SUTURE SIZE  
SUTURE NEEDLE  
MATERIAL




PRODUCT CODE  
SUTURE SIZE  
SUTURE NEEDLE  
MATERIAL



## SUTURE NEEDLES



### Securos 300™ Premium Suture Needles

Our needles are made from a propriety steel blend which is mechanically processed to make the material harder while still retaining high tensile strength. This enables the needles to have very sharp points and great durability while maintaining their flexibility. The same steel blend in our needles is used in the demanding human cardiac and plastic surgery areas where it has displaced more conventional needles.



### NEEDLE SHARPNESS & STRENGTH

Our needles are made from 300 series stainless steel, which has greater tensile strength, higher ductility and less brittleness than 400 series stainless steel.



Data on File

# NEEDLE TYPES

## REVERSE CUTTING ▽

3/8 Circle

C-3  
13mm 

FS-2  
19mm 

FS-1  
24mm 

FS  
26mm 

FSL  
30mm 

FSLX  
40mm 

LR  
70mm 

1/2 Circle

X-1  
22mm 

CP-2  
26mm 

CP-1  
36mm 

CP  
40mm 

## TAPERSHARP™ ▽

1/2 Circle

V-7  
26mm 

## PREMIUM REVERSE CUTTING ▽

3/8 Circle

P-1  
11mm 


P-3  
13mm 

PS-3  
16mm 

PS-2  
19mm 

## CONVENTIONAL CUTTING ▲

1/2 Circle

CCS-1  
60mm 

## STRAIGHT CUTTING ▽

KS  
60mm 

## SPATULA ▽

3/8 Circle

TG140-8  
6.5mm 

1/2 Circle

CS160-8  
5.5mm 

## TAPER POINT ⊙

1/2 Circle

TF  
13mm 

RB-1  
17mm 

SH-1  
22mm 


CT-3  
22mm 

SH  
26mm 


CT-2  
26mm 

CT-1  
36mm 

CTX  
48mm 



TP-1  
65mm 

## COBRA BLACK® NEEDLES



The Cobra Black® needle is a patent pending design that enhances visibility especially in blood and smaller surgery sites. The Cobra Black® needles eliminate the hindering glare from the operating room lighting. The needles are more distinguishable than the standard stainless steel needles.

### COBRA BLACK® vs. STANDARD NEEDLES



A hand wearing a white surgical glove holds a pair of surgical forceps. A blue suture thread is threaded through the forceps. The background is a plain, light-colored surface.

**ABSORBABLE SUTURE**  
**SECTION**

## NEEDLE TYPE

## SUTURE TYPE

Needle	Material	Color	Length	Suture Size (USP)					
				9-0	8-0	7-0	6-0	5-0	4-0
<b>XX-XX XXmm</b> X/X circle	XXXXX	XXXXX	<b>XX"</b> (XXcm)	SKU *					

♦ Cobra Black® Needle

## Premium Reverse Cutting ▽

## SECUROSORB™ QUICK (PGA)

Needle	Material	Color	Length	Suture Size					
				6-0	5-0	4-0	3-0	2-0	0
<b>P-1 11mm</b> 3/8 circle	Braided	Undyed	<b>18"</b> (45cm)	VGQ490 VGQ490C ♦					
<b>P-3 13mm</b> 3/8 circle	Braided	Undyed	<b>18"</b> (45cm)	VGQ494 VGQ494C					
<b>PS-3 16mm</b> 3/8 circle	Braided	Undyed	<b>30"</b> (75cm)	VGQ500C ♦ VGQ501C ♦					
<b>PS-2 19mm</b> 3/8 circle	Braided	Undyed	<b>18"</b> (45cm)	VGQ497 VGQ497C ♦					

♦ Cobra Black® Needle



## Taper Point ⊙

## SECUROSORB™ QUICK (PGA)

Needle	Material	Color	Length	Suture Size					
				6-0	5-0	4-0	3-0	2-0	0
<b>TF 13mm</b> 1/2 circle	Braided	Undyed	<b>18"</b> (45cm)	VGQ434C ♦					
<b>RB-1 17mm</b> 1/2 circle	Braided	Undyed	<b>18"</b> (45cm)	VGQ214C ♦					



♦ Cobra Black® Needle

**Spatula** 
**SECUROSORB™ (PGLA)**

					Suture Size					
Needle	Material	Color	Length	9-0	8-0	7-0	6-0	5-0	4-0	
<b>CS160-8 5.5mm</b> 3/8 circle	 Braided	Violet	<b>12"</b> (30cm)	<b>VG449*</b>						
<b>TG140-8 6.5mm</b> 3/8 circle	 Braided	Violet	<b>12"</b> (30cm)	<b>V548</b>						
		Violet	<b>18"</b> (45cm)	<b>V546</b>						







\*VG449 - Made with PGA suture material

**Premium Reverse Cutting** 
**SECUROSORB™ (PGLA)**

					Suture Size						
Needle	Material	Color	Length	7-0	6-0	5-0	4-0	3-0	2-0	0	
<b>P-1 11mm</b> 3/8 circle	 Braided	Undyed	<b>18"</b> (45cm)	<b>V488</b>	<b>V489</b>	<b>V490</b>					
<b>P-3 13mm</b> 3/8 circle	 Braided	Violet	<b>18"</b> (45cm)	<b>V463</b>							
		Undyed	<b>18"</b> (45cm)	<b>V494</b>							

### Reverse Cutting ▾


### SECUROSORB™ (PGLA)

				Suture Size						
Needle	Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2
<b>C-3 19mm</b> 3/8 circle		Violet	<b>18"</b> (45cm)	<b>V385</b>						
<b>FS-2 19mm</b> 3/8 circle		Braided	Violet	<b>18"</b> (45cm)	<b>V391</b>	<b>V392</b>				
			Violet	<b>30"</b> (75cm)		<b>V397</b> <b>V397C</b>	<b>V398</b> <b>V398C</b>			
<b>FS-1 24mm</b> 3/8 circle		Braided	Violet	<b>30"</b> (75cm)		<b>V451</b>	<b>V452</b>	<b>V453</b>		
			Undyed	<b>30"</b> (75cm)				<b>V452C</b> ♦ <b>V453C</b> ♦		
<b>FSL 30mm</b> 3/8 circle		Braided	Violet	<b>30"</b> (75cm)			<b>V586</b>	<b>V587</b>		
<b>X-1 22mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)			<b>V460</b>	<b>V461</b>		
<b>CP-1 36mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)			<b>V466</b>	<b>V467</b>	<b>V468</b>	
<b>CP 40mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)						<b>V195</b>

♦ Cobra Black® Needle

### Conventional Cutting ▲

### SECUROSORB™ (PGLA)

				Suture Size						
Needle	Material	Color	Length	4-0	3-0	2-0	0	1	2	3
<b>CCS-1 60mm</b> 1/2 circle	Braided	Violet	<b>30"</b> (75cm)							<b>V808</b>
										

**Taper Point**
**SECUROSORB™ (PGLA)**

					Suture Size						
Needle		Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2
<b>RB-1 17mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)	<b>V303</b>	<b>V304</b>					
<b>SH-1 22mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)		<b>V310</b>	<b>V311</b>				
<b>CT-3 22mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)			<b>V327</b>	<b>V328</b>			
<b>SH 26mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)			<b>V316</b>	<b>V317</b>			
							<b>V316C♦</b>				
<b>CT-2 26mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)			<b>V332</b>	<b>V333</b>	<b>V334</b>		
<b>CT-1 36mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)			<b>V338</b>	<b>V339</b>	<b>V340</b>	<b>V341</b>	
			Violet	<b>36"</b> (90cm)				<b>V345</b>	<b>V346</b>		
<b>CTX 48mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)						<b>V371</b>	
<b>TP-1 65mm</b> 1/2 circle		Braided	Violet	<b>54"</b> (137cm)							<b>V880</b>




♦ Cobra Black® Needle

**TaperSharp™**
**SECUROSORB™ (PGLA)**

					Suture Size						
Needle		Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2
<b>V-7 26mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)			<b>V9180</b>				





### Premium Reverse Cutting ▾

### SECUROCRYL™ (PGCL)

Needle	Material	Color	Length	Suture Size						
				6-0	5-0	4-0	3-0	2-0	0	1
<b>P-1 11mm</b> 3/8 circle 	Monofilament	Undyed	<b>18"</b> (45cm)	<b>Q489</b>						
<b>P-3 13mm</b> 3/8 circle 	Monofilament	Violet	<b>18"</b> (45cm)	<b>Q463</b>		<b>Q464</b>				
		Undyed	<b>18"</b> (45cm)	<b>Q493</b>		<b>Q494</b>				
<b>PS-2 19mm</b> 3/8 circle 	Monofilament	Undyed	<b>18"</b> (45cm)			<b>Q496</b>				





### Reverse Cutting ▾

### SECUROCRYL™ (PGCL)

Needle	Material	Color	Length	Suture Size						
				5-0	4-0	3-0	2-0	0	1	2
<b>FS-2 19mm</b> 3/8 circle 	Monofilament	Violet	<b>18"</b> (45cm)	<b>Q392</b>						
	Monofilament	Violet	<b>36"</b> (90cm)	<b>Q922</b>		<b>Q923</b>				
<b>FS-1 24mm</b> 3/8 circle 	Monofilament	Violet	<b>36"</b> (90cm)	<b>Q922C*</b>		<b>Q923C*</b>				
				<b>Q942</b>		<b>Q943</b>				
<b>FSL 30mm</b> 3/8 circle 	Monofilament	Violet	<b>36"</b> (90cm)					<b>Q987</b>		
				<b>Q942C*</b>		<b>Q943C*</b>				
<b>CP-1 36mm</b> 1/2 circle 	Monofilament	Violet	<b>36"</b> (90cm)					<b>Q967</b>		

\* Cobra Black® Needle


**Taper Point** 
**SECUROCRYL™ (PGCL)**

					Suture Size						
Needle		Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2
<b>RB-1 17mm</b> 1/2 circle		Monofilament	Violet	<b>30"</b> (75cm)	<b>Q303</b>	<b>Q304</b>					
					<b>Q303C♦</b>	<b>Q304C♦</b>					
		Monofilament	Undyed	<b>30"</b> (75cm)	<b>Q214</b>						
<b>SH 26mm</b> 1/2 circle		Monofilament	Violet	<b>30"</b> (75cm)	<b>Q315</b>		<b>Q316</b>				
								<b>Q317C♦</b>			
			Violet	<b>36"</b> (90cm)	<b>Q732</b>						
<b>CT-2 26mm</b> 1/2 circle		Monofilament	Violet	<b>30"</b> (75cm)	<b>Q334</b>						
								<b>Q762</b>			
			Violet	<b>36"</b> (90cm)							
<b>CT-1 36mm</b> 1/2 circle		Monofilament	Violet	<b>36"</b> (90cm)			<b>Q344</b>	<b>Q345</b>	<b>Q346</b>	<b>Q347</b>	

♦ Cobra Black® Needle






### Premium Reverse Cutting ▾

### SECURODOX™ (PDO)

Needle	Material	Color	Length	Suture Size						
				5-0	4-0	3-0	2-0	0	1	2
<b>P-3 13mm</b> 3/8 circle 	Monofilament	Violet	<b>18"</b> (45cm)	<b>D463</b>						







### Reverse Cutting ▾

### SECURODOX™ (PDO)

Needle	Material	Color	Length	Suture Size							
				5-0	4-0	3-0	2-0	0	1	2	
<b>FS-2 19mm</b> 3/8 circle 	Monofilament	Violet	<b>30"</b> (75cm)		<b>D397</b>	<b>D398</b>					
					<b>D397C♦</b>	<b>D398C♦</b>					
<b>FS-1 24mm</b> 3/8 circle 	Monofilament	Violet	<b>30"</b> (75cm)			<b>D452</b>	<b>D451</b>				
						<b>D452C♦</b>	<b>D451C♦</b>				
<b>CP-2 26mm</b> 1/2 circle 	Monofilament	Violet	<b>30"</b> (75cm)			<b>D968</b>	<b>D969</b>	<b>D970</b>			
<b>CP-1 36mm</b> 1/2 circle 	Monofilament	Violet	<b>30"</b> (75cm)				<b>D466</b>	<b>D467</b>	<b>D468</b>		
										<b>D467C♦</b>	
<b>CP 40mm</b> 1/2 circle 	Monofilament	Violet	<b>30"</b> (75cm)								<b>D195</b>

♦ Cobra Black® Needle


**Taper Point** 
**SECURODOX™ (PDO)**

Needle	Material	Color	Length	Suture Size						
				5-0	4-0	3-0	2-0	0	1	2
<b>RB-1</b> 17mm 1/2 circle 	Monofilament	Violet	<b>30"</b> (75cm)	<b>D303</b>	<b>D304</b>	<b>D305</b>				
<b>SH-1</b> 22mm 1/2 circle 	Monofilament	Violet	<b>30"</b> (75cm)		<b>D310</b>	<b>D311</b>				
<b>SH</b> 26mm 1/2 circle 	Monofilament	Violet	<b>30"</b> (75cm)		<b>D315</b>	<b>D316</b>	<b>D317</b> <b>D317C ♦</b>			
<b>CT-2</b> 26mm 1/2 circle 	Monofilament	Violet	<b>30"</b> (75cm)			<b>D332</b>	<b>D333</b>	<b>D334</b>		
<b>CT-1</b> 36mm 1/2 circle 	Monofilament	Violet	<b>30"</b> (75cm)			<b>D338</b>	<b>D339</b>	<b>D340</b>	<b>D341</b>	
<b>CT</b> 40mm 1/2 circle 	Monofilament	Violet	<b>36"</b> (90cm)					<b>D358</b>		

♦ Cobra Black® Needle


### Premium Reverse Cutting ▾

### CHROMIC GUT

Needle	Material	Color	Length	Suture Size						
				5-0	4-0	3-0	2-0	0	1	2
<b>P-3 13mm</b> 3/8 circle 	Twisted	Beige	<b>18"</b> (45cm)	CG1654						





### Reverse Cutting ▾

### CHROMIC GUT

Needle	Material	Color	Length	Suture Size							
				5-0	4-0	3-0	2-0	0	1	2	
<b>FS-2 19mm</b> 3/8 circle 	Twisted	Beige	<b>30"</b> (75cm)	CG634	CG635	CG636					
<b>FS-1 24mm</b> 3/8 circle 	Twisted	Beige	<b>30"</b> (75cm)				CG637				
<b>FS 26mm</b> 3/8 circle 	Twisted	Beige	<b>30"</b> (75cm)						CG638		
<b>CP-2 26mm</b> 1/2 circle 	Twisted	Beige	<b>30"</b> (75cm)				CG886	CG887			
<b>CP-1 36mm</b> 1/2 circle 	Twisted	Beige	<b>30"</b> (75cm)						CG817		

### Taper Point ⊙

### CHROMIC GUT

Needle	Material	Color	Length	Suture Size							
				5-0	4-0	3-0	2-0	0	1	2	
<b>SH-1 22mm</b> 1/2 circle 	Twisted	Beige	<b>30"</b> (75cm)	CG181							
<b>SH 26mm</b> 1/2 circle 	Twisted	Beige	<b>30"</b> (75cm)				CG122	CG123			
<b>CT-2 26mm</b> 1/2 circle 	Twisted	Beige	<b>30"</b> (75cm)				CG882	CG883	CG884		
<b>CT-1 36mm</b> 1/2 circle 	Twisted	Beige	<b>30"</b> (75cm)				CG811	CG812	CG813	CG814	

### Pre-Cut Sutures

### CHROMIC GUT

Needle	Material	Color	Length	Suture Size						
				3-0	2-0	0	1	2	3	4
<b>NONE</b>	Twisted	Beige	<b>54"</b> (135cm)	CG117						

# ANTIBACTERIAL SUTURE

PLUS SURGICAL SUTURE PRODUCTS

**A GREAT ALTERNATIVE FOR HIGH-RISK PATIENTS OR PROCEDURES WHERE BACTERIA IS A CONCERN.**

- ⊕ Securos PLUS sutures provide the same handling and performance of our other suture products, but with the added benefit of antibacterial protection.
- ⊕ Securos PLUS sutures are coated with chlorhexidine diacetate, a proven antibacterial agent, that result in a zone of inhibition around the suture line.
- ⊕ Our PLUS antibacterial agent is recognized as being effective against bacteria such as Staphylococcus aureus, Staphylococcus epidermidis, MRSA and MRSE.

## VETERINARY USE ONLY

### SECUROSORB™ PLUS



Coated, Braided, Polyglactin 910 (PGLA), Synthetic Absorbable Suture

Violet

### SECURODOX™ PLUS



Monofilament, Polydioxanone (PDO) Synthetic Absorbable Suture

Violet


### SECUROCRYL™ PLUS



Monofilament, Poly(glycolide-co-caprolactone (PGCL) Synthetic Absorbable Suture

Violet

## TENSILE STRENGTH CHART



## ZONE OF INHIBITION




PLUS Suture product with Antibacterial coating

\*Petri dish images is for illustrative purposes only, zone of inhibition testing results can vary



Suture product without Antibacterial coating

\*Petri dish images is for illustrative purposes only, zone of inhibition testing results can vary



MRSE - Vilet Plus Suture: size 2/0






S. aureus - Vilet Plus Suture: size 2/0

Images of Securos PLUS Suture product creating a zone of inhibition in vitro, reducing Staph aureus & MRSE bacterial colonization.



## Reverse Cutting ▾

## SECUROSORB™ PLUS (PGLA)

Antibacterial					Suture Size						
Needle		Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2
<b>FS-2 19mm</b> 3/8 circle		Braided	Undyed	<b>30"</b> (75cm)	<b>VP422</b>						
<b>FS-1 24mm</b> 3/8 circle		Braided	Undyed	<b>30"</b> (75cm)	<b>VP442</b>						
			Violet	<b>30"</b> (75cm)	<b>VP452</b>	<b>VP453</b>					
<b>CP-1 36mm</b> 1/2 circle		Braided	Undyed	<b>30"</b> (75cm)	<b>VP267</b>						
			Violet	<b>30"</b> (75cm)	<b>VP466</b>						




## Taper Point ⊙

## SECUROSORB™ PLUS (PGLA)

Antibacterial					Suture Size							
Needle		Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2	
<b>SH-1 22mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)	<b>VP310</b>							
<b>SH 26mm</b> 1/2 circle		Braided	Violet	<b>30"</b> (75cm)			<b>VP316</b>	<b>VP317</b>				




### Reverse Cutting ▾

### SECURODOX™ PLUS (PDO)

Antibacterial					Suture Size						
Needle		Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2
<b>FS-1 24mm</b> 3/8 circle		Monofilament	Violet	<b>30"</b> (75cm)			<b>DP452</b>	<b>DP451</b>			
<b>CP-2 26mm</b> 1/2 circle		Monofilament	Violet	<b>30"</b> (75cm)						<b>DP969</b>	
<b>CP-1 36mm</b> 1/2 circle		Monofilament	Violet	<b>30"</b> (75cm)						<b>DP466</b>	<b>DP467</b>




### Taper Point ⊙

### SECURODOX™ PLUS (PDO)

Antibacterial					Suture Size						
Needle		Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2
<b>SH 26mm</b> 1/2 circle		Monofilament	Violet	<b>30"</b> (75cm)			<b>DP315</b>	<b>DP316</b>			
<b>CT-2 26mm</b> 1/2 circle		Monofilament	Violet	<b>30"</b> (75cm)						<b>DP333</b>	
<b>CT-1 36mm</b> 1/2 circle		Monofilament	Violet	<b>30"</b> (75cm)						<b>DP339</b>	<b>DP340</b>



### Reverse Cutting ▾

### SECUROCRYL™ PLUS (PGCL)

Antibacterial					Suture Size					
Needle	Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2
<b>FS-2 19mm</b> 3/8 circle 	Monofilament	Violet	<b>36"</b> (90cm)		QP922	QP923				
<b>FS-1 24mm</b> 3/8 circle 	Monofilament	Violet	<b>36"</b> (90cm)			QP942	QP943			
<b>CP-1 36mm</b> 1/2 circle 	Monofilament	Violet	<b>36"</b> (90cm)							QP967

### Taper Point ⊙

### SECUROCRYL™ PLUS (PGCL)

Antibacterial					Suture Size					
Needle	Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2
<b>SH 26mm</b> 1/2 circle 	Monofilament	Violet	<b>30"</b> (75cm)			QP316	QP317			
<b>CT-1 36mm</b> 1/2 circle 	Monofilament	Violet	<b>36"</b> (90cm)			QP344	QP345	QP346		

**NON-ABSORBABLE SUTURE**  
**SECTION**




### Reverse Cutting ▾

### SECUROLON™ (NYLON)

Needle	Material	Color	Length	Suture Size								
				5-0	4-0	3-0	2-0	0	1	2		
<b>LR 75mm</b> 3/8 circle 	Monofilament	Blue	<b>20"</b> (50cm)									<b>460BL</b>

### Straight Cutting ▾

### SECUROLON™ (NYLON)

Needle	Material	Color	Length	Suture Size						
				7-0	5-0	4-0	3-0	2-0	0	
<b>KS 60mm</b> 	Monofilament	Black	<b>30"</b> (75cm)							<b>628BK</b>
		Blue	<b>30"</b> (75cm)							<b>628BL</b>


### Taper Point ⊙

### SECUROLON™ (NYLON)

Needle	Material	Color	Length	Suture Size							
				5-0	4-0	3-0	2-0	0	1	2	
<b>CTX 48mm</b> 1/2 circle 	Monofilament	Black	<b>30"</b> (75cm)								<b>830BK</b>
<b>TP-1 65mm</b> 1/2 circle 	Monofilament	Black	<b>60"</b> (150cm)								<b>825BK</b>

### Premium Reverse Cutting ▾






### SECUROLON™ (NYLON)

Needle	Material	Color	Length	Suture Size						
				7-0	6-0	5-0	4-0	3-0	2-0	0
<b>P-3 13mm</b> 3/8 circle 	Monofilament	Black	<b>18"</b> (45cm)	<b>698BK</b>						
		Blue								

◆ Cobra Black® Needle

### Reverse Cutting ▾


### SECUROLON™ (NYLON)

Needle	Material	Color	Length	Suture Size						
				6-0	5-0	4-0	3-0	2-0	0	1
<b>FS-2 19mm</b> 3/8 circle 	Monofilament	Black	<b>18"</b> (45cm)	<b>661BK</b>		<b>662BK</b>				
		Blue		<b>662BL</b>						
<b>FS-1 24mm</b> 3/8 circle 	Monofilament	Black	<b>18"</b> (45cm)	<b>663BK</b>						
		Blue		<b>663BL</b>						
<b>FS 26mm</b> 3/8 circle 	Monofilament	Black	<b>18"</b> (45cm)	<b>664BK</b>						
		Blue		<b>664BL</b>						
<b>FSL 30mm</b> 3/8 circle 	Monofilament	Blue	<b>30"</b> (75cm)	<b>665BL</b>						
<b>FSLX 40mm</b> 3/8 circle 	Monofilament	Blue	<b>30"</b> (75cm)	<b>1674BL</b>						

◆ Cobra Black® Needle






### Premium Reverse Cutting ▾

### SECUROPRO™ (POLYPROPYLENE)

Needle	Material	Color	Length	Suture Size							
				7-0	6-0	5-0	4-0	3-0	2-0	0	
<b>P-1 11mm</b> 3/8 circle 	Monofilament	Blue	<b>18"</b> (45cm)		<b>P8697</b> 8697						






### Reverse Cutting ▾

### SECUROPRO™ (POLYPROPYLENE)

Needle	Material	Color	Length	Suture Size							
				6-0	5-0	4-0	3-0	2-0	0	1	
<b>FS-2 19mm</b> 3/8 circle 	Monofilament	Blue	<b>18"</b> (45cm)		<b>P8661</b>	<b>P8683</b>	<b>P8665</b>				
<b>FS-1 24mm</b> 3/8 circle 	Monofilament	Blue	<b>18"</b> (45cm)				<b>P8684</b>				
<b>FS 26mm</b> 3/8 circle 	Monofilament	Blue	<b>18"</b> (45cm)					<b>P8685</b>			
<b>FSL 30mm</b> 3/8 circle 	Monofilament	Blue	<b>30"</b> (75cm)				<b>P8675</b>				
<b>FSLX 40mm</b> 3/8 circle 	Monofilament	Blue	<b>30"</b> (75cm)					<b>P8689</b>	<b>P8690</b>		





## Taper Point

## SECUROPRO™ (POLYPROPYLENE)

					Suture Size						
Needle		Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2
<b>RB-1 17mm</b> 1/2 circle		Monofilament	Blue	<b>30"</b> (75cm)		<b>P8871</b>	<b>P8872</b>				
<b>SH 26mm</b> 1/2 circle		Monofilament	Blue	<b>30"</b> (75cm)			<b>P8832</b>				
<b>CT-2 26mm</b> 1/2 circle		Monofilament	Blue	<b>30"</b> (75cm)				<b>P8411</b>			
<b>CT-1 36mm</b> 1/2 circle		Monofilament	Blue	<b>30"</b> (75cm)				<b>P8423</b>	<b>P8424</b>	<b>P8425</b>	
<b>TP-1 65mm</b> 1/2 circle		Monofilament	Blue	<b>60"</b> (150cm)							<b>P8825</b>

## Reverse Cutting ▾

## VISIPRO\*\* (POLYPROPYLENE)

Needle	Material	Color	Length	Suture Size						
				5-0	4-0	3-0	2-0	0	1	
<b>FS-2 19mm</b> 3/8 circle 	Monofilament	Fluorescent Pink	<b>30"</b> (75cm)		<b>F8683</b>	<b>F8665</b>				
<b>FS-1 24mm</b> 3/8 circle 	Monofilament	Fluorescent Pink	<b>30"</b> (75cm)			<b>F8684</b>				
<b>FS 26mm</b> 3/8 circle 	Monofilament	Fluorescent Pink	<b>30"</b> (75cm)				<b>F8685</b>			
<b>FSLX 24mm</b> 3/8 circle 	Monofilament	Fluorescent Pink	<b>30"</b> (75cm)				<b>F8689</b>	<b>F8690</b>		

\*\* VISIPRO is a trademark of Riverpoint Medical

**Premium Reverse Cutting ▾**
**SECUROSILK™ (SILK)**

					Suture Size						
Needle		Material	Color	Length	7-0	6-0	5-0	4-0	3-0	2-0	0
<b>P-1 11mm</b> 3/8 circle		Braided	Black	<b>18"</b> (45cm)	<b>S639</b>						
<b>P-3 13mm</b> 3/8 circle		Braided	Black	<b>18"</b> (45cm)	<b>S640</b>						

**Reverse Cutting ▾**
**SECUROSILK™ (SILK)**

					Suture Size						
Needle		Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2
<b>C-3 13mm</b> 3/8 circle		Braided	Black	<b>18"</b> (45cm)	<b>S736</b>						
<b>FS-2 19mm</b> 3/8 circle		Braided	Black	<b>18"</b> (45cm)	<b>S683</b>						
<b>FS-1 24mm</b> 3/8 circle		Braided	Black	<b>18"</b> (45cm)	<b>S684</b>						
<b>FS 26mm</b> 3/8 circle		Braided	Black	<b>18"</b> (45cm)	<b>S685</b>						
<b>FSL 30mm</b> 3/8 circle		Braided	Black	<b>30"</b> (75cm)	<b>S680</b>						

**Taper Point ⊕**
**SECUROSILK™ (SILK)**

					Suture Size						
Needle		Material	Color	Length	5-0	4-0	3-0	2-0	0	1	2
<b>SH 26mm</b> 1/2 circle		Braided	Black	<b>30"</b> (75cm)	<b>S833</b>						

# SUTURE CASSETTES SECTION



# SUTURE MICRO-CASSETTES

## CASSETTE FEATURES

EASY-OPEN LID

CONVENIENT BUILT-IN CUTTER

ERGONOMICAL DESIGN

POCKET SIZE

ECO-FRIENDLY  
(50% less plastic than traditional cassette packaging)



## CASSETTE RACK



SKU	DESCRIPTION
MCR-100	Cassette Rack

Fits 5 Suture Cassettes

SECUROS SECUROSORB™	SUTURE MATERIAL	SUTURE COLOR	SUTURE SIZE (metric)	SUTURE LENGTH
MCV-3	Braided	Violet	3 (6)	25 meters
MCV-2/0	Braided	Violet	2/0 (3.0)	25 meters
MCV-3/0	Braided	Violet	3/0 (2.0)	25 meters

SECUROS SECURODOX™	SUTURE MATERIAL	SUTURE COLOR	SUTURE SIZE (metric)	SUTURE LENGTH
MCD-2	Monofilament	Violet	2 (5)	25 meters
MCD-1	Monofilament	Violet	1 (4)	25 meters
MCD-0	Monofilament	Violet	0 (3.5)	25 meters
MCD-2/0	Monofilament	Violet	2/0 (3.0)	25 meters
MCD-3/0	Monofilament	Violet	3/0 (2.0)	25 meters
MCD-4/0	Monofilament	Violet	4/0 (1.5)	25 meters

SECUROS SECUROCRYL™	SUTURE MATERIAL	SUTURE COLOR	SUTURE SIZE (metric)	SUTURE LENGTH
MCQ-0	Monofilament	Violet	0 (3.5)	25 meters
MCQ-2/0	Monofilament	Violet	2/0 (3.0)	25 meters
MCQ-3/0	Monofilament	Violet	3/0 (2.0)	25 meters

SECUROS VISIPRO**	SUTURE MATERIAL	SUTURE COLOR	SUTURE SIZE (metric)	SUTURE LENGTH
MCF-0	Monofilament	Pink	0 (3.5)	85 meters
MCF-2/0	Monofilament	Pink	2/0 (3.0)	110 meters
MCF-3/0	Monofilament	Pink	3/0 (2.0)	110 meters
MCF-4/0	Monofilament	Pink	4/0 (1.5)	110 meters

SECUROS CHROMIC GUT	SUTURE MATERIAL	SUTURE COLOR	SUTURE SIZE (metric)	SUTURE LENGTH
MCG-4	Twisted	Beige	4 (8)	25 meters
MCG-3	Twisted	Beige	3 (7)	40 meters
MCG-2	Twisted	Beige	2 (6)	45 meters
MCG-1	Twisted	Beige	1 (5)	50 meters
MCG-0	Twisted	Beige	0 (4)	50 meters
MCG-2/0	Twisted	Beige	2/0 (3.5)	50 meters
MCG-3/0	Twisted	Beige	3/0 (3)	50 meters

SECUROS SECURUMID™	SUTURE MATERIAL	SUTURE COLOR	SUTURE SIZE (metric)	SUTURE LENGTH
MCS-4	Braided	White	4 (6)	25 meters
MCS-3	Braided	White	3 (6)	50 meters
MCS-2	Braided	White	2 (5)	50 meters
MCS-1	Braided	White	1 (4)	75 meters
MCS-0	Braided	White	0 (3.5)	75 meters
MCS-2/0	Braided	White	2/0 (3.0)	100 meters
MCS-3/0	Braided	White	3/0 (2.0)	100 meters
MCS-4/0	Braided	White	4/0 (1.5)	100 meters

\*\* VISIPRO is a trademark of Riverpoint Medical



## MATERIAL COMPARISON CHARTS

### ABSORBABLE SUTURE

RIVERPOINT	ETHICON	COVIDIEN	DAVIS & GECK	VPL
<b>SECUROSORB™ QUICK</b>	Vicryl Rapide <sup>^</sup>	•	•	•
<b>SECUROSORB™</b>	Vicryl <sup>^</sup>	Polysorb <sup>^^</sup>	Dexon <sup>^^^</sup>	Polymend ES <sup>^^^^</sup>
<b>SECUROCRYL™</b>	Monocryl <sup>^</sup>	Caprosyn/Biosyn <sup>^^</sup>	•	Monomend MT <sup>^^^^</sup>
<b>SECURODOX™</b>	PDS II <sup>^</sup>	•	•	Monomend Max <sup>^^^^</sup>
<b>CHROMIC GUT</b>	Chromic Gut	Chromic Gut	Chromic Gut	Chromed <sup>^^^^</sup>
<b>SECUROSORB™ PLUS</b>	Vicryl Plus <sup>^</sup>	•	•	•
<b>SECUROCRYL™ PLUS</b>	Monocryl Plus <sup>^</sup>	•	•	•
<b>SECURODOX™ PLUS</b>	PDS Plus <sup>^</sup>	•	•	•

### NON-ABSORBABLE SUTURE

RIVERPOINT	ETHICON	COVIDIEN	DAVIS & GECK	VPL
<b>SECUROLON™</b>	Ethilon <sup>^</sup>	Monosof <sup>^^</sup>	Dermalon <sup>^^^</sup>	NY-STAN <sup>^^^^</sup>
<b>SECUROPRO™</b>	Prolene <sup>^</sup>	Surgipro <sup>^^</sup>	•	PRO-STAN <sup>^^^^</sup>
<b>SECUROSILK™</b>	Silk	Silk	Silk	•
<b>VISIPRO™</b>	•	•	•	•

<sup>^</sup>trademark of Ethicon Inc.

<sup>^^</sup>trademark of Covidien

<sup>^^^</sup>trademark of Surgical Specialties Corp.

<sup>^^^^</sup>trademark of Riverpoint Medical

# SKU QUICK REFERENCE CHARTS—ORDERING

## ABSORBABLE SUTURE

SECUROSORB™				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
V808	039256	3	CCS-1	30"
V195	039259	2	CP	30"
V880	039265	2	TP-1	54"
V341	039241	1	CT-1	30"
V371	065658	1	CTX	30"
V468	038173	1	CP-1	30"
V587	039225	0	FSL	30"
V334	039235	0	CT-2	30"
V340	038178	0	CT-1	30"
V346	065668	0	CT-1	36"
V467	038172	0	CP-1	30"
V443	065665	2-0	FS-1	30"
V453	038180	2-0	FS-1	30"
V453C ♦	062716	2-0	FS-1	30"
V586	039255	2-0	FSL	30"
V461	039266	2-0	X-1	30"
V328	065676	2-0	CT-3	30"
V317	038179	2-0	SH	30"
V333	038177	2-0	CT-2	30"
V339	039224	2-0	CT-1	30"
V345	039260	2-0	CT-1	36"
V466	038171	2-0	CP-1	30"
V398	038181	3-0	FS-2	30"
V398C ♦	062719	3-0	FS-2	30"
V452	038175	3-0	FS-1	30"
V452C ♦	062717	3-0	FS-1	30"
V460	039257	3-0	X-1	30"
V311	039261	3-0	SH-1	30"
V327	039267	3-0	CT-3	30"
V316	039231	3-0	SH	30"
V316C ♦	062720	3-0	SH	30"
V332	038176	3-0	CT-2	30"
V338	065732	3-0	CT-1	30"
V9180	370243	3-0	V-7	30"
V392	065753	4-0	FS-2	18"
V397	038174	4-0	FS-2	30"
V397C ♦	062718	4-0	FS-2	30"
V451	039238	4-0	FS-1	30"
V494	065759	4-0	P-3	18"
V304	065793	4-0	RB-1	30"
V310	039227	4-0	SH-1	30"
V385	065711	5-0	C-3	18"
V391	065726	5-0	FS-2	18"
V490	039275	5-0	P-1	18"
V463	039272	5-0	P-3	18"
V303	065765	5-0	RB-1	30"
V489	039258	6-0	P-1	18"

SECUROSORB™				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
V488	065702	7-0	P-1	18"
V546	981334	7-0	TG140-8	18"
V548	039277	8-0	TG140-8	12"
VG449	056012	9-0	CS160-8	12"

SECUROSORB™ QUICK				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
VGQ497	059274	3-0	PS-2	18"
VGQ497C ♦	062698	3-0	PS-2	18"
VGQ494	059275	4-0	P-3	18"
VGQ494C ♦	062699	4-0	P-3	18"
VGQ501C ♦	062703	4-0	PS-3	30"
VGQ434C ♦	062704	4-0	TF	18"
VGQ214C ♦	062701	4-0	RB-1	18"
VGQ490	059276	5-0	P-1	18"
VGQ490C ♦	062700	5-0	P-1	18"
VGQ500C ♦	062702	5-0	PS-3	30"

SECUROCRYL™				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
Q347	038160	1	CT-1	36"
Q987	038167	0	FSL	36"
Q334	065758	0	CT-2	30"
Q346	038158	0	CT-1	36"
Q967	038155	0	CP-1	36"
Q943	038163	2-0	FS-1	36"
Q943C ♦	062705	2-0	FS-1	36"
Q732	038162	2-0	SH	36"
Q317C ♦	062709	2-0	SH	30"
Q762	038159	2-0	CT-2	36"
Q345	038166	2-0	CT-1	36"
Q923	038164	3-0	FS-2	36"
Q923C ♦	062707	3-0	FS-2	36"
Q942	038157	3-0	FS-1	36"
Q942C ♦	062706	3-0	FS-1	36"
Q316	038161	3-0	SH	30"
Q344	039233	3-0	CT-1	36"
Q392	065169	4-0	FS-2	18"
Q922	038156	4-0	FS-2	36"
Q922C ♦	062708	4-0	FS-2	36"
Q464	039273	4-0	P-3	18"
Q494	062678	4-0	P-3	18"
Q496	065756	4-0	PS-2	18"
Q214	065791	4-0	RB-1	30"
Q304	064382	4-0	RB-1	30"
Q304C ♦	065843	4-0	RB-1	30"

SECUROCRYL™				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
Q315	038165	4-0	SH	30"
Q463	039239	5-0	P-3	18"
Q493	065773	5-0	P-3	18"
Q303	064381	5-0	RB-1	30"
Q303C ♦	064383	5-0	RB-1	30"
Q489	065686	6-0	P-1	18"

SECUROCRYL™				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
D195	039254	2	CP	30"
D341	038145	1	CT-1	30"
D468	039234	1	CP-1	30"
D334	039228	0	CT-2	30"
D340	038144	0	CT-1	30"
D467	039148	0	CP-1	30"
D467C ♦	062715	0	CP-1	30"
D970	038138	0	CP-2	30"
D358	065755	0	CT	36"
D451	038135	2-0	FS-1	30"
D451C ♦	062713	2-0	FS-1	30"
D317	038147	2-0	SH	30"
D317C ♦	062710	2-0	SH	30"
D333	038146	2-0	CT-2	30"
D339	038143	2-0	CT-1	30"
D466	038137	2-0	CP-1	30"
D969	039222	2-0	CP-2	30"
D398	038072	3-0	FS-2	30"
D398C ♦	062712	3-0	FS-2	30"
D452	038073	3-0	FS-1	30"
D452C ♦	062714	3-0	FS-1	30"
D305	039245	3-0	RB-1	30"
D311	038140	3-0	SH-1	30"
D316	038141	3-0	SH	30"
D332	039221	3-0	CT-2	30"
D338	038142	3-0	CT-1	30"
D968	038136	3-0	CP-2	30"
D397	038071	4-0	FS-2	30"
D397C ♦	062711	4-0	FS-2	30"
D304	039237	4-0	RB-1	30"
D310	038139	4-0	SH-1	30"
D315	038989	4-0	SH	30"
D463	049119	5-0	P-3	18"
D303	065772	5-0	RB-1	30"

♦ Cobra Black® Needle

# SKU QUICK REFERENCE CHARTS—ORDERING

## ABSORBABLE SUTURE

CHROMIC GUT				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
CG117	065684	3	NONE	54"
CG814	065687	2	CT-1	30"
CG813	039262	1	CT-1	30"
CG638	065692	0	FS	30"
CG884	039268	0	CT-2	30"
CG812	038152	0	CT-1	30"
CG817	065683	0	CP-1	30"
CG637	038151	2-0	FS-1	30"
CG123	038149	2-0	SH	30"
CG883	038154	2-0	CT-2	30"
CG811	065698	2-0	CT-1	30"
CG887	065688	2-0	CP-2	30"
CG636	038150	3-0	FS-2	30"
CG122	038148	3-0	SH	30"
CG882	038153	3-0	CT-2	30"
CG886	065653	3-0	CP-2	30"
CG635	039252	4-0	FS-2	30"
CG1654	065675	4-0	P-3	18"
CG181	039263	4-0	SH-1	30"
CG634	065731	5-0	FS-2	30"

SECUROSORB™ PLUS				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
VP267	065790	0	CP-1	30"
VP453	064903	2-0	FS-1	30"
VP317	065787	2-0	SH	30"
VP466	065657	2-0	CP-1	30"
VP442	065701	3-0	FS-1	30"
VP452	064902	3-0	FS-1	30"
VP316	065727	3-0	SH	30"
VP422	065789	4-0	FS-2	30"
VP310	065656	4-0	SH-1	30"

ANTIBACTERIAL

SECUROCRYL™ PLUS				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
QP346	065767	0	CT-1	36"
QP967	065792	0	CP-1	36"
QP943	064910	2-0	FS-1	36"
QP317	065771	2-0	SH	30"
QP345	064906	2-0	CT-1	36"
QP923	064907	3-0	FS-2	36"
QP942	064908	3-0	FS-1	36"
QP316	064909	3-0	SH	30"
QP344	065754	3-0	CT-1	36"
QP922	065770	4-0	FS-2	36"

ANTIBACTERIAL

SECURODOX™ PLUS				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
DP340	065690	0	CT-1	30"
DP467	065785	0	CP-1	30"
DP451	064899	2-0	FS-1	30"
DP333	065786	2-0	CT-2	30"
DP339	065689	2-0	CT-1	30"
DP466	065764	2-0	CP-1	30"
DP969	065708	2-0	CP-2	30"
DP452	064905	3-0	FS-1	30"
DP316	064904	3-0	SH	30"
DP315	065763	4-0	SH	30"

ANTIBACTERIAL

## NON-ABSORBABLE SUTURE

SECUROLON™				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
460BL	039269	2	LR	20"
825BK	065685	2	TP-1	60"
830BK	065700	1	CTX	30"
665BL	049117	0	FSL	30"
664BK	048359	2-0	FS	18"
664BL	038170	2-0	FS	30"
1674BL	039226	2-0	FSLX	30"
628BK	065762	2-0	KS	30"
628BL	049118	2-0	KS	30"
663BK	090060	3-0	FS-1	18"
663BL	038169	3-0	FS-1	30"
662BK	048358	4-0	FS-2	18"
662BL	038168	4-0	FS-2	30"
661BK	045210	5-0	FS-2	18"
698BK	065729	5-0	P3	18"

SECUROPRO™				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
P8825	039274	2	TP-1	60"
P8425	065757	1	CT-1	30"
P8690	039253	0	FSLX	30"
P8424	039264	0	CT-1	30"
P8685	039230	2-0	FS	18"
P8689	065709	2-0	FSLX	30"
P8411	039271	2-0	CT-2	30"
P8423	065654	2-0	CT-1	30"
P8665	039244	3-0	FS-2	18"
P8684	039243	3-0	FS-1	18"
P8675	065672	3-0	FSL	30"
P8872	065678	3-0	RB-1	30"
P8832	042124	3-0	SH	30"
P8683	039242	4-0	FS-2	18"
P8871	065677	4-0	RB-1	30"
P8661	065695	5-0	FS-2	18"
P8697	045064	6-0	P-1	18"

SECUROSILK™				
SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
S680	065681	0	FSL	30"
S685	980929	2-0	FS	18"
S833	065682	2-0	SH	30"
S684	039276	3-0	FS-1	18"
S736	065691	4-0	C-3	18"
S683	039270	4-0	FS-2	18"
S640	065693	5-0	P-3	18"
S639	065699	6-0	P-1	18"

VISI PRO™

SECUROS ITEM CODE	MWI ITEM CODE	SIZE	NEEDLE	LENGTH
F8690	039285	0	FSLX	30"
F8685	039284	2-0	FS	30"
F8689	039286	2-0	FSLX	30"
F8665	039287	3-0	FS-2	30"
F8684	039283	3-0	FS-1	30"
F8683	039282	4-0	FS-2	30"

\*\* VISI PRO is a trademark of Riverpoint Medical

# SKU QUICK REFERENCE CHARTS—COMPARISON

## ABSORBABLE SUTURE

SECUROS SECUROSORB™	ETHICON Vicryl <sup>^</sup>
V195	J195H
V303	J303H
V304	J304H
V310	J310H
V311	J311H
V316	J316H
V316C ♦	J316H
V317	J317H
V327	J327H
V328	J328H
V332	J332H
V333	J333H
V334	J334H
V338	J338H
V339	J339H
V340	J340H
V341	J341H
V345	J345H
V346	J346H
V371	J371H
V385	J385H
V391	J391H
V392	J392H
V397	J397H
V397C ♦	J397H
V398	J398H
V398C ♦	J398H
V443	J443H
VG449 ♦	J449H

SECUROS SECUROSORB™	ETHICON Vicryl <sup>^</sup>
V451	J451H
V452	J452H
V452C ♦	J452H
V453	J453H
V453C ♦	J453H
V460	J460H
V461	J461H
V463	J463G
V466	J466H
V467	J467H
V468	J468H
V488	J488H
V489	J489G
V490	J490G
V494	J494G
V546	J546H
V548	J548H
V586	J586H
V587	J587H
V808	J808T
V880	J880T
V9180	W9180

SECUROS SECURODOX™	ETHICON PDS III <sup>^</sup>
D195	Z195H
D303	Z303H
D304	Z304H
D305	Z305H
D310	Z310H
D311	Z311H
D315	Z315H
D316	Z316H
D317	Z317H
D317C ♦	Z317H
D332	Z332H
D333	Z333H
D334	Z334H
D338	Z338H
D339	Z339H
D340	Z340H
D341	Z341H
D358	Z358H
D397	Z397H
D397C ♦	Z397H
D398	Z398H
D398C ♦	Z397H
D451	Z451H
D451C ♦	Z397H
D452	Z452H
D452C ♦	Z452H
D463	Z463H
D466	Z466H
D467	Z467H
D467C ♦	Z467H
D468	Z468H
D968	Z968H
D969	Z969H
D970	Z970H

SECUROS SECUROCRYL™	ETHICON Monocryl <sup>^</sup>
Q214	Y214H
Q303	Y303H
Q303C ♦	•
Q304	Y304H
Q304C ♦	•
Q310	Y310H
Q315	Y315H
Q316	Y316H
Q317C ♦	•
Q334	Y334H
Q344	Y344H
Q345	Y345H
Q346	Y346H
Q347	Y347H
Q392	Y392H
Q463	Y463G
Q464	Y464G
Q489	Y489G
Q493	Y493G
Q494	Y494G
Q496	Y496G
Q732	Y732H
Q762	Y762H
Q922	Y922H
Q922C ♦	•
Q923	Y923H
Q923C ♦	•
Q942	Y942H
Q942C ♦	Y942H
Q943	Y943H

SECUROS SECUROCRYL™	ETHICON Monocryl <sup>^</sup>
Q943C ♦	•
Q967	Y967H
Q987	Y987H

SECUROS SECUROSORB™ QUICK	ETHICON Vicryl Rapide <sup>^</sup>
VGQ214C ♦	VR214
VGQ434C ♦	VR434
VGQ490	VR490
VGQ490C ♦	•
VGQ494	VR494
VGQ494C ♦	•
VGQ497	VR497
VGQ497C ♦	•
VGQ500C ♦	•
VGQ501C ♦	•

♦ Cobra Black® Needle

<sup>^</sup>trademark of Ethicon Inc.

# SKU QUICK REFERENCE CHARTS—COMPARISON

## ABSORBABLE SUTURE

SECUROS CHROMIC GUT	ETHICON Chromic Gut
CG117	S117H
CG122	G122H
CG123	G123H
CG181	G181H
CG634	634G
CG635	635G
CG636	636G
CG637	637G
CG638	638G
CG811	811H
CG812	812H
CG813	813H
CG814	814H
CG817	817H
CG882	882H
CG883	883H
CG884	884H
CG886	886H
CG887	887H
CG1654	1654H

SECUROS SECUROSORB™ PLUS	ETHICON Vicryl Plus <sup>^</sup>
VP267	VCP267H
VP310	VCP310H
VP316	VCP316H
VP317	VCP317H
VP422	VCP422H
VP442	VCP442H
VP452	VCP52H
VP453	VCP453H
VP466	VCP466H
<b>ANTIBACTERIAL</b>	

SECUROS SECURODOX™ PLUS	ETHICON PDS Plus <sup>^</sup>
DP315	PDP315H
DP316	PDP316H
DP333	PDP333H
DP339	PDP339H
DP340	PDP340H
DP451	PDP451H
DP452	PDP452H
DP466	PDP466H
DP467	PDP467H
DP969	PDP452H
<b>ANTIBACTERIAL</b>	

SECUROS SECUROCRYL™ PLUS	ETHICON Monocryl Plus <sup>^</sup>
QP316	MCP316H
QP317	MCP317H
QP344	MCP344H
QP345	MCP345H
QP346	MCP346H
QP922	MCP922H
QP923	MCP923H
QP942	MCP942H
QP943	MCP943H
QP967	MCP967H
<b>ANTIBACTERIAL</b>	

<sup>^</sup>trademark of Ethicon Inc.

## NON-ABSORBABLE SUTURE

SECUROS SECUROLON™	ETHICON Ethilon <sup>^</sup>
460BL	460H
628BK	628H
628BL	628H
661BK	661G
662BK	662G
662BL	662G
663BK	663G
663BL	663G
664BK	664G
664BL	664G
665BL	665H
698BK	698G
825BK	825H
830BK	830H
1674BL	1674H

SECUROS SECUROPRO™	ETHICON Prolene <sup>^</sup>
P8411	8411H
P8423	8423H
P8424	8424H
P8425	8425H
P8661	8661G
P8665	8665G
P8675	8675H
P8683	8683G
P8684	8684G
P8685	8685H
P8687	8687H
P8689	8689H
P8690	8690H
P8697	8697G
P8825	8825H

SECUROS SECUROPRO™	ETHICON Prolene <sup>^</sup>
P8697	8697G
P8825	8825H
P8832	8832H
P8871	8871H
P8872	8872H

SECUROS SECUROSILK™	ETHICON Silk
S639	639G
S640	640G
S680	680G
S683	683G
S684	684G
S685	685G
S736	736G
S833	833G

SECUROS VISIPRO**	ETHICON •
F8665	•
F8683	•
F8684	•
F8685	•
F8689	•
F8690	•

<sup>^</sup>trademark of Ethicon Inc.  
<sup>\*\*</sup>trademark of Riverpoint Medical



Where knowledge,  
reach and partnership  
shape healthcare delivery.