

SawCover System

The SawCover System consists of a cover and a hardware saw that is optimized with an ergonomic handle. Paired together, they allow you to saw with ease and confidence.

The ANSI Level 4 medical textile cover meets all required standards of EN13795, ASTM1670, ASTM1671 and therefore is impermeable to blood, bodily fluids, bacteria and water.

Peace of Mind

Feel confident knowing that your saw is completely sterile thanks to the highest level medical barrier on the market.

Keep Orthopedics In-House

The SawCover and DrillCover Systems together are the affordable power tools to get started in orthopedics.

Item #	Description	Arbutus Catalog #
066100	SAW KIT COMPLETE ARBUTUS	KVI-002 (SW-A101, SYS-005)
066101	SAW COVER SYSTEM ARBUTUS	SYS -005 (RP-009, RP-010, CHK-003)
066102	Sagittal Saw Blade Arbutus 20mm	BLD - 201
066103	Sagittal Saw Blade Arbutus 13mm	BLD - 202
066104	Sagittal Saw Blade Arbutus 10mm	BLD - 203
066106	Sag Saw Blade Variety Pk 2x BLD-203, 2x BLD-202, 1x BLD-201	BLD - V200



Inner fabric:
Woven polyester layer protects ePTFE membrane.

ePTFE membrane:
Prevents transfer of liquids, bacteria, pathogens, viruses.

Outer fabric: Woven polyester layer for added protection and durability.





Ergonomic pistol grip handle offers easy handling of saw in any position with extra support and leverage



Securely tighten your blade with the thumb screw



Continuously adjustable speed from 5,000 to 20,000 RPM



Instant ON/OFF trigger switch for improved patient and operator safety



Saw Blades

A range of blades are available for all procedures involving an oscillating saw: osteotomy, THA, TKA, trauma and reconstruction, amputation.

Blades are made of the highest quality stainless steel by the same facilities in the USA and Europe that make blades for major orthopaedic device companies.



Linen Specifications

Barrier Level	Level 4 ANSI/AAMI PB70 barrier performance against fluid penetration
Lifetime	75 use cycles
Registration	Arbutus Medical is a registered establishment with the U.S. FDA and is ISO 13485 certified
Sterilization	Steam autoclave 30 mins exposure time at 121°C

Saw Specifications

Voltage	12 V battery, rechargeable on either 110-120 V or 220-240 V AC (charger dependent)
Speed	Variable speed (5,000 - 20,000 OPM)
Battery	2x lithium ion, rechargeable, 30 minutes continuous unloaded runtime on full charge
Ergonomic	Weighs only 1.1kg and is easy to handle due to a custom fitted pistol-grip handle
Warranty	1 year manufacturer's warranty

Arbutus Medical

#560-828 West 10th Ave
Vancouver, BC V5Z1L8



Arbutus Customer Service

 veterinary@arbutusmedical.ca

 +1-778-296-2577 Option #4