

Securos Surgical®

Stryker®

Power equipment



Trusted, reliable power at affordable prices

Securos now offers preowned, Stryker brand, human-grade drilling and sawing surgical power equipment as an economical alternative to new, full-price equipment. Stryker equipment was developed by working closely with orthopedic surgeons to deliver solutions to the complex challenges of today's operating rooms. If you'd like a backup set of equipment for your practice or something reliable and lower-cost for everyday use, Stryker is what you're looking for in trusted, reliable power.

The equipment features:

- Leading ergonomics to give surgeons neutral wrist positioning, allowing for firm and positive grip
- Components, including batteries, are autoclaveable for greater efficiency between procedures
- Equipment kits available in complete drill and saw packages and individual drill or saw only options
- Lightweight, rechargeable batteries that power the handpieces
- 2 or 4 port battery charging systems
- 2 new Stryker or Economy battery options

Handpiece kit options



Stryker 1 Handpiece, Drill Kit

123717

Includes one of each (preowned)

- CD4 Cordless Driver/Handpiece
- Stryker Pin Collet
- Stryker Wire Collet
- Stryker AO Attachment
- Stryker 1/4" Jacobs Chuck and Key
- Stryker 5/32" Jacobs Chuck and Key



Stryker 2 Handpiece, Drill and Saw Kit

123718

Includes one of each (preowned)

- CD4 Cordless Driver/Handpiece
- Sabo2 Sagittal Saw Driver/Handpiece
- Stryker Pin Collet
- Stryker Wire Collet
- Stryker AO Attachment
- Stryker 1/4" Jacobs Chuck and Key
- Stryker 5/32" Jacobs Chuck and Key

Warranty and service:

Stryker being a longtime leader in orthopedic power equipment, we are confident that you will love your Stryker system. Hardware items come with a one year warranty and batteries with a six month warranty along with Securos's trusted and reliable customer service and support.

Confidence in your power

The cordless CD4 driver handpiece gives you confidence for drilling and cutting needs during stifle remediation and orthopedic trauma procedures.


**This kit includes these features
and quick-change attachments:**

- Stryker dual-trigger, CD4 handpiece for forward, reverse, and oscillating rotation
- Separate pin and wire collets
- 1/4" keyed drill chuck
- 5/32" keyed drill chuck
- AO Quick Connect chuck

Your saw cutting requirements are addressed with either a separate Stryker Sabo2 Saw handpiece or an optional Sagittal Saw attachment. Either package is lightweight and delivers optimized cutting speed, control, and accuracy during surgery.



Charging system options



Stryker Charging Kit

123849

- 1 Stryker 4-port charging station (preowned)
- 2 small, Li-Ion batteries (new)



Economy Charging Kit

123848

- 1 economy 2-port charging station (new)
- 2 large NiCad batteries (new)

Sagittal saw blade options



Stryker Sagittal Saw Attachment

123713

- Sold separately or included in 123844 - Handpiece kit

Stryker Sagittal Saw Blades

123972	13 x 39 x 0.38mm
124005	9 x 25x 0.38mm
123968	9 x 31 x 0.38mm

Europe +49 7467 9476 50
UK +44 845 070 2498
USA +1 877.266.3349 option #2

info@securos-europe.eu
sales@securos.co.uk
securossales@mwiah.com

Premium **quality**,
uncompromising **value**

securos.com