

SurgX

Interlocking Nail

Instrument and implant kits
composition catalog




Securos Surgical®

INSTRUMENT KITS


CATS & SMALL DOGS

NAILS 3.5 / 4.0 / 5.0 SHORT

Cats & small dogs instrument tray



Common instrument tray #1




Common instrument tray #2



Sterilization Boxes




Common instrument tray n°1



Rep.	Description	Qty
1	Handle with AO quick coupling	1
2	Square awl reamer for nails \varnothing 3.5 / 4 / 5	1
3	Guide wire \varnothing 1.2 for modular intramedullary reamer size 1	1
4	Flexible shaft for modular intramedullary reamer size 1	1
5	Reamer head Size 1 - \varnothing 3.5	1
	Reamer head Size 1 - \varnothing 4	1
	Reamer head Size 1 - \varnothing 5	1
	Reamer head Size 1 - \varnothing 6	1
6	Depth gauge Size 1 for screws \varnothing 2.0 - 2.5	1
7	Depth gauge Size 2 for screws \varnothing 3.0 - 3.5	1
8	T8 star driver for screws \varnothing 2.0-2.5	1
9	T15 star driver for screws \varnothing 3.0-3.5	1
10	Soft tissue protector - Size 1*	1
11	Common instruments sterilization tray #1	1


* Optional instrument not included in basic set

Common instrument tray n°2



Rep.	Description	Qty
1	Hex screwdriver 3.5	1
2	Bracket for targeting system	1
3	Bracket wheel-Anvil	1
4	Targeting arm wheel	1
5	Drill Sleeve	4
6	CIS-cortex Drill Bit Ø2.4	1
7	TRANS-cortex Drill Bit Ø1.5 Size 1	2
8	TRANS-cortex Drill Bit Ø1.5 Size 2	2
9	CIS-cortex Drill Bit Ø3.0	1
10	TRANS-cortex Drill Bit Ø2.0 Size 1	2
11	TRANS-cortex Drill Bit Ø2.0 Size 2	2
12	CIS-cortex Drill Bit Ø4.0	1
13	TRANS-cortex Drill Bit Ø2.4 Size 1	2
14	TRANS-cortex Drill Bit Ø2.4 Size 2	2
15	Common instruments sterilization tray #2	1
-	Common instruments sterilization container - 310 x 140 x 105mm - Grey lid	1

Cats & Small Dogs Instrument Tray



Rep.	Description	Qty
1	Nail extension Size 1 Ø3.5	1
2	Nail extension Size 1 Ø4.0	1
3	Nail extension Size 1 Ø5.0	1
4	Nail extension screw Size 1	1
5	Guide wire Ø1.2 Lg 270mm	1
6	Targeting Arm Size 1 for nails Ø3.5 / 4 / 5	1
7	Cats & small dogs instrument sterilization tray - Nails 3.5 - 4 - 5	1
-	Cats & small dogs instrument sterilization container - 310 x 140 x 76mm - Yellow lid	1


MEDIUM & LARGE DOGS

NAILS 5.0 LONG / 6.0 / 8.0

Medium & large dogs instrument tray



Common instrument tray #1




Common instrument tray #2



Sterilization Boxes




Common instrument tray n°1



Rep.	Description	Qty
1	Handle with AO quick coupling	1
2	Square awl reamer for nails Ø5 / 6 / 8	1
3	Guide wire Ø1.5 for modular intramedullary reamer size 2	1
4	Flexible shaft for modular intramedullary reamer size 2	1
5	Reamer head Size 2 - Ø5	1
	Reamer head Size 2 - Ø6	1
	Reamer head Size 2 - Ø7	1
	Reamer head Size 2 - Ø8	1
6	Reamer head Size 2 - Ø9	1
	Depth gauge Size 2 for screws Ø3.0 - 3.5	1
7	T15 star driver for screws Ø3.0-3.5	1
8	Soft tissue protector - Size 2 *	1
9	Common instruments sterilization tray #1	1

* Optional instrument not included in basic set

Common instrument tray n°2




Rep.	Description	Qty
1	Hex screwdriver 3.5	1
2	Bracket for targeting system	1
3	Bracket wheel-Anvil	1
4	Targeting arm wheel	1
5	Drill Sleeve	4
6	CIS-cortex Drill Bit Ø4.0	1
7	TRANS-cortex Drill Bit Ø2.4 Size 1	2
8	TRANS-cortex Drill Bit Ø2.4 Size 2	2
9	TRANS-cortex Drill Bit Ø2.8 Size 1	2
10	TRANS-cortex Drill Bit Ø2.8 Size 2	2
11	Common instruments sterilization tray #2	1
-	Common instruments sterilization container - 310 x 140 x 105mm - Grey lid	1

CATS & DOGS – ALL SIZES

NAILS 3.5 / 4.0 / 5.0 SHORT / 5.0 LONG / 6.0 / 8.0


Cats & small dogs instrument tray



Medium & large dogs instrument tray



Common instrument tray #1




Common instrument tray #2



Sterilization Boxes




Medium & Large Dogs Instrument Tray



Rep.	Description	Qty
1	Nail extension Size 2 Ø5.0	1
2	Nail extension Size 2 Ø6.0	1
3	Nail extension Size 2 Ø8.0	1
4	Nail extension screw Size 2	1
5	Nitinol guide wire Ø1.5 Lg 350mm	1
6	Targeting Arm Size 2 - Proximal part	1
7	Targeting Arm Size 2 - Distal part	1
8	Targeting Arm Size 2 - Wheel	1
9	Medium & large dogs instrument sterilization tray - Nails 5 - 8	1
-	Medium & large dogs instrument sterilization container - 310 x 140 x 76mm - Blue lid	1


Common instrument tray n°1



Rep.	Description	Qty
1	Handle with AO quick coupling	1
2	Square awl reamer for nails Ø3.5 / 4 / 5	1
3	Square awl reamer for nails Ø5 / 6 / 8	1
4	Guide wire Ø1.2 for modular intramedullary reamer size 1	1
5	Guide wire Ø1.5 for modular intramedullary reamer size 2	1
6	Flexible shaft for modular intramedullary reamer size 1	1
7	Flexible shaft for modular intramedullary reamer size 2	1
8	Reamer head Size 1 - Ø3.5	1
	Reamer head Size 1 - Ø4	1
	Reamer head Size 1 - Ø5	1
	Reamer head Size 1 - Ø6	1
9	Reamer head Size 2 - Ø5	1
	Reamer head Size 2 - Ø6	1
	Reamer head Size 2 - Ø7	1
	Reamer head Size 2 - Ø8	1
10	Depth gauge Size 1 for screws Ø2.0 - 2.5	1
	Depth gauge Size 2 for screws Ø3.0 - 3.5	1
	T8 star driver for screws Ø2.0-2.5	1
12	T8 star driver for screws Ø2.0-2.5	1
13	T15 star driver for screws Ø3.0-3.5	1
14	Soft tissue protector - Size 1 *	1
15	Soft tissue protector - Size 2 *	1
16	Common instruments sterilization tray #1	1


* Optional instruments not included in basic set

Common instrument tray n°2




Rep.	Description	Qty
1	Hex screwdriver 3.5	1
2	Bracket for targeting system	1
3	Bracket wheel-Anvil	1
4	Targeting arm wheel	1
5	Drill Sleeve	4
6	CIS-cortex Drill Bit Ø2.4	1
7	TRANS-cortex Drill Bit Ø1.5 Size 1	2
8	TRANS-cortex Drill Bit Ø1.5 Size 2	2
9	CIS-cortex Drill Bit Ø3.0	1
10	TRANS-cortex Drill Bit Ø2.0 Size 1	2
11	TRANS-cortex Drill Bit Ø2.0 Size 2	2
12	CIS-cortex Drill Bit Ø4.0	1
13	TRANS-cortex Drill Bit Ø2.4 Size 1	2
14	TRANS-cortex Drill Bit Ø2.4 Size 2	2
15	TRANS-cortex Drill Bit Ø2.8 Size 1	2
16	TRANS-cortex Drill Bit Ø2.8 Size 2	2
17	Common instruments sterilization tray #2	1
-	Common instruments sterilization container - 310 x 140 x 105mm - Grey lid	1

Cats & Small Dogs Instrument Tray



Rep.	Description	Qty
1	Nail extension Size 1 Ø3.5	1
2	Nail extension Size 1 Ø4.0	1
3	Nail extension Size 1 Ø5.0	1
4	Nail extension screw Size 1	1
5	Guide wire Ø1.2 Lg 270mm	1
6	Targeting Arm Size 1 for nails Ø3.5 / 4 / 5	1
7	Cats & small dogs instrument sterilization tray - Nails 3.5 - 4 - 5	1
-	Cats & small dogs instrument sterilization container - 310 x 140 x 76mm - Yellow lid	1

Medium & Large Dogs Instrument Tray




Rep.2	Description	Qty
1	Nail extension Size 2 Ø5.0	1
2	Nail extension Size 2 Ø6.0	1
3	Nail extension Size 2 Ø8.0	1
4	Nail extension screw Size 2	1
5	Nitinol guide wire Ø1.5 Lg 350mm	1
6	Targeting Arm Size 2 - Proximal part	1
7	Targeting Arm Size 2 - Distal part	1
8	Targeting Arm Size 2 - Wheel	1
9	Medium & large dogs instrument sterilization tray - Nails 5 - 6 - 8	1
-	Medium & large dogs instrument sterilization container - 310 x 140 x 76mm - Blue lid	1

IMPLANT KITS


CATS & SMALL DOGS

NAILS 3.5 / 4.0 / 5.0 SHORT

Nails 3.5 / 4.0 tray



Nails 5.0 s tray



Screw rack




Sterilization Boxes




IMPLANT SET COMPOSITION
INTERLOCKING NAIL SYSTEM
CATS & SMALL DOGS - STARTER SET

Interlocking Nails ϕ 3.5 / 4.0




Ref.	Description	Qty
1	Interlocking Nail ϕ 3.5 Lg 70mm	1
	Interlocking Nail ϕ 3.5 Lg 78mm	1
	Interlocking Nail ϕ 3.5 Lg 86mm	1
	Interlocking Nail ϕ 3.5 Lg 94mm	1
	Interlocking Nail ϕ 3.5 Lg 102mm	1
	Interlocking Nail ϕ 3.5 Lg 110mm	1
2	Interlocking Nail ϕ 4.0 Lg 86mm	1
	Interlocking Nail ϕ 4.0 Lg 94mm	1
	Interlocking Nail ϕ 4.0 Lg 102mm	1
3	Interlocking Nail ϕ 4.0 Lg 110mm	1
-	Sterilization tray for nails ϕ 3.5 / 4.0	1
-	Sterilization container for nails ϕ 3.5 / 4 - 310 x 140 x 45mm - Yellow lid	1

Interlocking Nails ϕ 5.0 short



Ref.	Description	Qty
	Interlocking Nail ϕ 5.0 Lg 94mm	1
	Interlocking Nail ϕ 5.0 Lg 102mm	1
	Interlocking Nail ϕ 5.0 Lg 110mm	1
4	Interlocking Nail ϕ 5.0 Lg 118mm	1
	Interlocking Nail ϕ 5.0 Lg 126mm	1
	Interlocking Nail ϕ 5.0 Lg 134mm	1
	Interlocking Nail ϕ 5.0 Lg 142mm	1
	Interlocking Nail ϕ 5.0 Lg 150mm	1
5	Sterilization tray for nails ϕ 5 short	1
-	Sterilization container for nails ϕ 5.0 - 310 x 140 x 45mm - Green lid	1

Screws & Bolts



Ref.	Description	Qty
	Screw for Interlocking Nail Ø2.0 Lg 10mm	3
	Screw for Interlocking Nail Ø2.0 Lg 12mm	3
6	Screw for Interlocking Nail Ø2.0 Lg 14mm	3
	Screw for Interlocking Nail Ø2.0 Lg 16mm	2
	Screw for Interlocking Nail Ø2.0 Lg 18mm	2
	Screw for Interlocking Nail Ø2.0 Lg 20mm	2
	Screw for Interlocking Nail Ø2.5 Lg 12mm	3
	Screw for Interlocking Nail Ø2.5 Lg 14mm	3
7	Screw for Interlocking Nail Ø2.5 Lg 16mm	3
	Screw for Interlocking Nail Ø2.5 Lg 18mm	3
	Screw for Interlocking Nail Ø2.5 Lg 20mm	2
	Screw for Interlocking Nail Ø2.5 Lg 22mm	2
	Screw for Interlocking Nail Ø2.5 Lg 24mm	2
	Screw for Interlocking Nail Ø2.5 Lg 26mm	2
	Screw for Interlocking Nail Ø3.0 Lg 14mm	2
	Screw for Interlocking Nail Ø3.0 Lg 16mm	3
	Screw for Interlocking Nail Ø3.0 Lg 18mm	3
	Screw for Interlocking Nail Ø3.0 Lg 20mm	3
	Screw for Interlocking Nail Ø3.0 Lg 22mm	3
8	Screw for Interlocking Nail Ø3.0 Lg 24mm	2
	Screw for Interlocking Nail Ø3.0 Lg 26mm	2
	Screw for Interlocking Nail Ø3.0 Lg 28mm	2
	Screw for Interlocking Nail Ø3.0 Lg 30mm	2
	Screw for Interlocking Nail Ø3.0 Lg 32mm	2
	Screw for Interlocking Nail Ø3.0 Lg 34mm	2
	Screw for Interlocking Nail Ø3.0 Lg 36mm	2


Screws & Bolts (Continued)

Ref.	Description	Qty
	Bolt for Interlocking Nail Ø2.0 Lg 10mm	1
	Bolt for Interlocking Nail Ø2.0 Lg 12mm	1
9	Bolt for Interlocking Nail Ø2.0 Lg 14mm	1
	Bolt for Interlocking Nail Ø2.0 Lg 16mm	1
	Bolt for Interlocking Nail Ø2.0 Lg 18mm	1
	Bolt for Interlocking Nail Ø2.0 Lg 20mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 12mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 14mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 16mm	1
10	Bolt for Interlocking Nail Ø2.5 Lg 18mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 20mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 22mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 24mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 26mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 14mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 16mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 18mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 20mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 22mm	1
11	Bolt for Interlocking Nail Ø3.0 Lg 24mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 26mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 28mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 30mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 32mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 34mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 36mm	1
12	Screw rack	1
-	Sterilization container for screw rack - 310 x 140 x 105mm - Grey lid	1

MEDIUM & LARGE DOGS

NAILS 5.0 LONG / 6.0 / 8.0


Nails 5.0 L tray



Nails 6.0 tray



Nails 8.0 tray




Screw rack



Sterilization Boxes




Interlocking Nails ϕ 5.0 long




Ref.	Description	Qty
1	Interlocking Nail ϕ 5.0 Lg 134mm	1
	Interlocking Nail ϕ 5.0 Lg 146.5mm	1
	Interlocking Nail ϕ 5.0 Lg 159mm	1
	Interlocking Nail ϕ 5.0 Lg 171.5mm	1
2	Sterilization tray for nails ϕ 5 long	1
-	Sterilization container for nails ϕ 5.0 - 310 x 140 x 45mm - Green lid	1

Interlocking Nails ø6.0




Ref.	Description	Qty
	Interlocking Nail Ø6.0 Lg 134mm	1
	Interlocking Nail Ø6.0 Lg 146.5mm	1
3	Interlocking Nail Ø6.0 Lg 159mm	1
	Interlocking Nail Ø6.0 Lg 171.5mm	1
	Interlocking Nail Ø6.0 Lg 184mm	1
4	Sterilization tray for nails ø6	1
-	Sterilization container for nails ø6.0 - 310 x 140 x 45mm - Blue lid	1

Interlocking Nails ø8.0



Ref.	Description	Qty
	Interlocking Nail Ø8.0 Lg 159mm	1
	Interlocking Nail Ø8.0 Lg 171.5mm	1
5	Interlocking Nail Ø8.0 Lg 184mm	1
	Interlocking Nail Ø8.0 Lg 196.5mm	1
6	Sterilization tray for nails ø8	1
-	Sterilization container for nails ø8.0 - 310 x 140 x 45mm - Purple lid	1

Screws & Bolts



Ref.	Description	Qty	
7	Screw for Interlocking Nail Ø3.0 Lg 14mm	2	
	Screw for Interlocking Nail Ø3.0 Lg 16mm	3	
	Screw for Interlocking Nail Ø3.0 Lg 18mm	3	
	Screw for Interlocking Nail Ø3.0 Lg 20mm	3	
	Screw for Interlocking Nail Ø3.0 Lg 22mm	3	
	Screw for Interlocking Nail Ø3.0 Lg 24mm	2	
	Screw for Interlocking Nail Ø3.0 Lg 26mm	2	
	Screw for Interlocking Nail Ø3.0 Lg 28mm	2	
	Screw for Interlocking Nail Ø3.0 Lg 30mm	2	
	Screw for Interlocking Nail Ø3.0 Lg 32mm	2	
	Screw for Interlocking Nail Ø3.0 Lg 34mm	2	
	Screw for Interlocking Nail Ø3.0 Lg 36mm	2	
	8	Screw for Interlocking Nail Ø3.5 Lg 16mm	2
		Screw for Interlocking Nail Ø3.5 Lg 18mm	2
Screw for Interlocking Nail Ø3.5 Lg 20mm		3	
Screw for Interlocking Nail Ø3.5 Lg 22mm		3	
Screw for Interlocking Nail Ø3.5 Lg 24mm		3	
Screw for Interlocking Nail Ø3.5 Lg 26mm		3	
Screw for Interlocking Nail Ø3.5 Lg 28mm		2	
Screw for Interlocking Nail Ø3.5 Lg 30mm		2	
Screw for Interlocking Nail Ø3.5 Lg 32mm		2	
Screw for Interlocking Nail Ø3.5 Lg 34mm		2	
Screw for Interlocking Nail Ø3.5 Lg 36mm		2	
Screw for Interlocking Nail Ø3.5 Lg 38mm		2	
Screw for Interlocking Nail Ø3.5 Lg 40mm		2	
Screw for Interlocking Nail Ø3.5 Lg 42mm		2	
Screw for Interlocking Nail Ø3.5 Lg 45mm	2		
Screw for Interlocking Nail Ø3.5 Lg 48mm	2		

IMPLANT SET COMPOSITION
INTERLOCKING NAIL SYSTEM
MEDIUM & LARGE DOGS - STARTER SET

REF

IS000210


Screws & Bolts (Continued)

Ref.	Description	Qty	
9	Bolt for Interlocking Nail Ø3.0 Lg 14mm	1	
	Bolt for Interlocking Nail Ø3.0 Lg 16mm	1	
	Bolt for Interlocking Nail Ø3.0 Lg 18mm	1	
	Bolt for Interlocking Nail Ø3.0 Lg 20mm	1	
	Bolt for Interlocking Nail Ø3.0 Lg 22mm	1	
	Bolt for Interlocking Nail Ø3.0 Lg 24mm	1	
	Bolt for Interlocking Nail Ø3.0 Lg 26mm	1	
	Bolt for Interlocking Nail Ø3.0 Lg 28mm	1	
	Bolt for Interlocking Nail Ø3.0 Lg 30mm	1	
	Bolt for Interlocking Nail Ø3.0 Lg 32mm	1	
	Bolt for Interlocking Nail Ø3.0 Lg 34mm	1	
	Bolt for Interlocking Nail Ø3.0 Lg 36mm	1	
	10	Bolt for Interlocking Nail Ø3.5 Lg 16mm	1
		Bolt for Interlocking Nail Ø3.5 Lg 18mm	1
Bolt for Interlocking Nail Ø3.5 Lg 20mm		1	
Bolt for Interlocking Nail Ø3.5 Lg 22mm		1	
Bolt for Interlocking Nail Ø3.5 Lg 24mm		1	
Bolt for Interlocking Nail Ø3.5 Lg 26mm		1	
Bolt for Interlocking Nail Ø3.5 Lg 28mm		1	
Bolt for Interlocking Nail Ø3.5 Lg 30mm		1	
Bolt for Interlocking Nail Ø3.5 Lg 32mm		1	
Bolt for Interlocking Nail Ø3.5 Lg 34mm		1	
Bolt for Interlocking Nail Ø3.5 Lg 36mm		1	
Bolt for Interlocking Nail Ø3.5 Lg 38mm		1	
Bolt for Interlocking Nail Ø3.5 Lg 40mm		1	
Bolt for Interlocking Nail Ø3.5 Lg 42mm		1	
Bolt for Interlocking Nail Ø3.5 Lg 45mm	1		
Bolt for Interlocking Nail Ø3.5 Lg 48mm	1		
11	Screw rack	1	
-	Sterilization container for screw rack - 310 x 140 x 105mm - Grey lid	1	


CATS & DOGS – ALL SIZES

NAILS 3.5 / 4.0 / 5.0 SHORT / 5.0 LONG / 6.0 / 8.0

Nails 3.5 / 4.0 tray



Nails 5.0 s tray




Nails 5.0 L tray



Nails 6.0 tray



Nails 8.0 tray




Screw rack



Sterilization Boxes




Interlocking Nails ϕ 3.5 / 4.0




Ref.	Description	Qty
1	Interlocking Nail ϕ 3.5 Lg 70mm	1
	Interlocking Nail ϕ 3.5 Lg 78mm	1
	Interlocking Nail ϕ 3.5 Lg 86mm	1
	Interlocking Nail ϕ 3.5 Lg 94mm	1
	Interlocking Nail ϕ 3.5 Lg 102mm	1
	Interlocking Nail ϕ 3.5 Lg 110mm	1
2	Interlocking Nail ϕ 4.0 Lg 86mm	1
	Interlocking Nail ϕ 4.0 Lg 94mm	1
	Interlocking Nail ϕ 4.0 Lg 102mm	1
	Interlocking Nail ϕ 4.0 Lg 110mm	1
3	Sterilization tray for nails ϕ 3.5 / 4.0	1
-	Sterilization container for nails ϕ 3.5 / 4 - 310 x 140 x 45mm - Yellow lid	1

Interlocking Nails ø5.0 short




Ref.	Description	Qty
4	Interlocking Nail Ø5.0 Lg 94mm	1
	Interlocking Nail Ø5.0 Lg 102mm	1
	Interlocking Nail Ø5.0 Lg 110mm	1
	Interlocking Nail Ø5.0 Lg 118mm	1
	Interlocking Nail Ø5.0 Lg 126mm	1
	Interlocking Nail Ø5.0 Lg 134mm	1
	Interlocking Nail Ø5.0 Lg 142mm	1
	Interlocking Nail Ø5.0 Lg 150mm	1
5	Sterilization tray for nails ø5 short	1
-	Sterilization container for nails ø5.0 - 310 x 140 x 45mm - Green lid	1

Interlocking Nails \varnothing 5.0 long




Ref.	Description	Qty
6	Interlocking Nail \varnothing 5.0 Lg 159mm	1
	Interlocking Nail \varnothing 5.0 Lg 171.5mm	1
7	Sterilization tray for nails \varnothing 5 long	1
-	Sterilization container for nails \varnothing 5.0 - 310 x 140 x 45mm - Green lid	1

Interlocking Nails ø6.0




Ref.	Description	Qty
	Interlocking Nail Ø6.0 Lg 134mm	1
	Interlocking Nail Ø6.0 Lg 146.5mm	1
8	Interlocking Nail Ø6.0 Lg 159mm	1
	Interlocking Nail Ø6.0 Lg 171.5mm	1
	Interlocking Nail Ø6.0 Lg 184mm	1
9	Sterilization tray for nails ø6	1
-	Sterilization container for nails ø6.0 - 310 x 140 x 45mm - Blue lid	1

Interlocking Nails ø8.0



Ref.	Description	Qty
	Interlocking Nail Ø8.0 Lg 159mm	1
10	Interlocking Nail Ø8.0 Lg 171.5mm	1
	Interlocking Nail Ø8.0 Lg 184mm	1
	Interlocking Nail Ø8.0 Lg 196.5mm	1
11	Sterilization tray for nails ø8	1
-	Sterilization container for nails ø8.0 - 310 x 140 x 45mm - Purple lid	1

Screws & Bolts



Ref.	Description	Qty	
12	Screw for Interlocking Nail Ø2.0 Lg 10mm	3	
	Screw for Interlocking Nail Ø2.0 Lg 12mm	3	
	Screw for Interlocking Nail Ø2.0 Lg 14mm	3	
	Screw for Interlocking Nail Ø2.0 Lg 16mm	2	
	Screw for Interlocking Nail Ø2.0 Lg 18mm	2	
	Screw for Interlocking Nail Ø2.0 Lg 20mm	2	
13	Screw for Interlocking Nail Ø2.5 Lg 12mm	3	
	Screw for Interlocking Nail Ø2.5 Lg 14mm	3	
	Screw for Interlocking Nail Ø2.5 Lg 16mm	3	
	Screw for Interlocking Nail Ø2.5 Lg 18mm	3	
	Screw for Interlocking Nail Ø2.5 Lg 20mm	2	
	Screw for Interlocking Nail Ø2.5 Lg 22mm	2	
	Screw for Interlocking Nail Ø2.5 Lg 24mm	2	
	Screw for Interlocking Nail Ø2.5 Lg 26mm	2	
	14	Screw for Interlocking Nail Ø3.0 Lg 14mm	2
		Screw for Interlocking Nail Ø3.0 Lg 16mm	3
Screw for Interlocking Nail Ø3.0 Lg 18mm		3	
Screw for Interlocking Nail Ø3.0 Lg 20mm		3	
Screw for Interlocking Nail Ø3.0 Lg 22mm		3	
Screw for Interlocking Nail Ø3.0 Lg 24mm		2	
Screw for Interlocking Nail Ø3.0 Lg 26mm		2	
Screw for Interlocking Nail Ø3.0 Lg 28mm		2	
Screw for Interlocking Nail Ø3.0 Lg 30mm		2	
Screw for Interlocking Nail Ø3.0 Lg 32mm		2	
Screw for Interlocking Nail Ø3.0 Lg 34mm		2	
Screw for Interlocking Nail Ø3.0 Lg 36mm		2	
15		Screw for Interlocking Nail Ø3.5 Lg 16mm	2
		Screw for Interlocking Nail Ø3.5 Lg 18mm	2
	Screw for Interlocking Nail Ø3.5 Lg 20mm	3	
	Screw for Interlocking Nail Ø3.5 Lg 22mm	3	
	Screw for Interlocking Nail Ø3.5 Lg 24mm	3	
	Screw for Interlocking Nail Ø3.5 Lg 26mm	3	
	Screw for Interlocking Nail Ø3.5 Lg 28mm	2	
	Screw for Interlocking Nail Ø3.5 Lg 30mm	2	
	Screw for Interlocking Nail Ø3.5 Lg 32mm	2	
	Screw for Interlocking Nail Ø3.5 Lg 34mm	2	
	Screw for Interlocking Nail Ø3.5 Lg 36mm	2	
	Screw for Interlocking Nail Ø3.5 Lg 38mm	2	
	Screw for Interlocking Nail Ø3.5 Lg 40mm	2	
	Screw for Interlocking Nail Ø3.5 Lg 42mm	2	
Screw for Interlocking Nail Ø3.5 Lg 45mm	2		
Screw for Interlocking Nail Ø3.5 Lg 48mm	2		

Screws & Bolts (Continued)

Ref.	Description	Qty
16	Bolt for Interlocking Nail Ø2.0 Lg 10mm	1
	Bolt for Interlocking Nail Ø2.0 Lg 12mm	1
	Bolt for Interlocking Nail Ø2.0 Lg 14mm	1
	Bolt for Interlocking Nail Ø2.0 Lg 16mm	1
	Bolt for Interlocking Nail Ø2.0 Lg 18mm	1
	Bolt for Interlocking Nail Ø2.0 Lg 20mm	1
17	Bolt for Interlocking Nail Ø2.5 Lg 12mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 14mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 16mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 18mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 20mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 22mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 24mm	1
	Bolt for Interlocking Nail Ø2.5 Lg 26mm	1
18	Bolt for Interlocking Nail Ø3.0 Lg 14mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 16mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 18mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 20mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 22mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 24mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 26mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 28mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 30mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 32mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 34mm	1
	Bolt for Interlocking Nail Ø3.0 Lg 36mm	1
19	Bolt for Interlocking Nail Ø3.5 Lg 16mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 18mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 20mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 22mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 24mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 26mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 28mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 30mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 32mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 34mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 36mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 38mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 40mm	1
	Bolt for Interlocking Nail Ø3.5 Lg 42mm	1
Bolt for Interlocking Nail Ø3.5 Lg 45mm	1	
Bolt for Interlocking Nail Ø3.5 Lg 48mm	1	
20	Screw rack	1
-	Sterilization container for screw rack - 310 x 140 x 105mm - Grey lid	1



Ø3.5

70 - 142mm



< 7 kg

Ø4.0

78 - 142mm



< 10 kg

Ø5.0

94 - 196.5mm



< 22 kg

Ø6.0

134 - 221.5mm



< 30 kg

Ø8.0

146.5 - 234mm



> 30 kg

Securos Surgical
Securos USA
433 Main Street
Fiskdale, MA 01518
securos.com

Tel: 877.266.3349
Email: info@securos.com

Securos Surgical
Securos Europe GmbH
Take off Gewerbepark 11
78579 Neuhausen ob Eck - Germany
securos.com

Tél: +49 07467 9489120
Email: info@securos-europe.eu

Surg'X
26 rue Barthélémy de Laffemas
Technosite
26000 Valence - France

Tél: +33 (0)4 75 25 88 21 Email:
contact@surgx.fr

Premium **quality**,
uncompromising **value**

[securos.com](https://www.secuross.com)