

Premium Suture

Premium Antibacterial "PLUS" Suture



Save Big


Securos Surgical can save you up to 76% on suture when switching from the industry leader. The average savings per hospital based on the most common suture offerings is 41%.

Ask your MWI Animal Health sales representative about our 100% money-back guarantee and sample the Securos Suture line today!

- A great alternative for high-risk patients or procedures where bacteria is a concern.
- Securos PLUS sutures provide the same handling and performance of our other suture products, but with the added benefit of antibacterial protection.
- Securos PLUS sutures are coated with chlorhexidine diacetate, a proven antibacterial agent, resulting in a zone of inhibition around the suture line.
- Our PLUS antibacterial agent is recognized as being effective against bacteria such as Staphylococcus aureus, Staphylococcus epidermidis, MRSA and MRSE.

Additional Features

- Racetrack packaging design delivers low suture memory
- One dozen packaging
- Proudly made in USA



SKUS ON NEXT PAGE



We build *our* work around the way *you* work.

Securos Surgical sold exclusively through MWI Animal Health

☎ 800.824.3703

☎ Service Support: 877.266.3349

🌐 mwiah.com

🌐 securos.com

Premium Antibacterial "PLUS" Suture

| CODE | NEEDLE | SIZE | LENGTH |
|-------------------------|--------|------|--------|
| SECUROCRYL™ Plus | | | |
| QP967 | CP-1 | 0 | 36" |
| QP346 | CT-1 | 0 | 36" |
| QP345 | CT-1 | 2-0 | 36" |
| QP943 | FS-1 | 2-0 | 36" |
| QP317 | SH | 2-0 | 30" |
| QP344 | CT-1 | 3-0 | 36" |
| QP942 | FS-1 | 3-0 | 36" |
| QP923 | FS-2 | 3-0 | 36" |
| QP316 | SH | 3-0 | 30" |
| QP922 | FS-2 | 4-0 | 36" |

| CODE | NEEDLE | SIZE | LENGTH |
|------------------------|--------|------|--------|
| SECURODOX™ Plus | | | |
| DP467 | CP-1 | 0 | 30" |
| DP340 | CT-1 | 0 | 30" |
| DP466 | CP-1 | 2-0 | 30" |
| DP969 | CP-2 | 2-0 | 30" |
| DP339 | CT-1 | 2-0 | 30" |
| DP333 | CT-2 | 2-0 | 30" |
| DP451 | FS-1 | 2-0 | 30" |
| DP452 | FS-1 | 3-0 | 30" |
| DP316 | SH | 3-0 | 30" |
| DP315 | SH | 4-0 | 30" |

| CODE | NEEDLE | SIZE | LENGTH |
|-------------------------|--------|------|--------|
| SECUROSORB™ Plus | | | |
| VP466 | CP-1 | 2-0 | 30" |
| VP453 | FS-1 | 2-0 | 30" |
| VP317 | SH | 2-0 | 30" |
| VP452 | FS-1 | 3-0 | 30" |
| VP442 | FS-1 | 3-0 | 30" |
| VP316 | SH | 3-0 | 30" |
| VP422 | FS-2 | 4-0 | 30" |
| VP310 | SH-1 | 4-0 | 30" |