

Premium Suture

VISIPRO Non-Absorbable Fluorescent Suture



Save Big

Securos Surgical can save you up to 74% on suture when switching from the industry leader. The average savings per hospital based on the most common suture offerings is 40%.


Ask your MWI Animal Health sales representative about our 100% money-back guarantee and sample the Securos Suture line today!

- Fluoresces under a black light
- Bright magenta color – great for dark coated animals
- Monofilament polypropylene for smooth tissue passage and excellent knot security

Additional Features

- Racetrack packaging design delivers low suture memory
- One dozen packaging
- Proudly made in USA

CODE	NEEDLE	SIZE	LENGTH
VISIPRO			
F8690	FSLX	0	30"
F8685	FS	2-0	30"
F8689	FSLX	2-0	30"
F8684	FS-1	3-0	30"
F8665	FS-2	3-0	30"
F8683	FS-2	4-0	30"
MCF-0	CASSETTE	0	85M
MCF-2/0	CASSETTE	2-0	110M
MCF-3/0	CASSETTE	3-0	110M
MCF-4/0	CASSETTE	4-0	110M



We build *our* work around the way *you* work.

Securos Surgical sold exclusively through MWI Animal Health

800.824.3703

Service Support: 877.266.3349

mwiah.com

securos.com